

Supplemental Table S1. Primers designed for qRT-PCR analysis

Gene	Forward primers (5'-3')	Reverse primers (5'-3')
<i>MsG0180001827.01</i>	AACCTGGATTTTCGCGGTGAT	TACACCCAACCTCCTCCCCT
<i>MsG0180005519.01</i>	TAGGGCAGTGGTAAACAGCG	CAGGGTAGATCCTTGGGTTGT
<i>MsG0180005358.01</i>	GAAAGCCACCAGAGAAGTGC	CCCACTCAAGTCCTTCTCCC
<i>MsG0780040428.01</i>	CCGGTAAAATGGCCGAAAGC	GAGCGGGGTTAGAGTGGATG
<i>MsG0780036812.01</i>	TGCCATCATCGGTAACAACCA	TGGGACGTGCGAAAAGAGTT
<i>MsG0480020076.01</i>	CAAGGTATGGTTTGCACCGC	GAGATCCTCCGATACCCCCA
<i>MsG0480018569.01</i>	CTGTGACCAGCACCATTCCT	CTTGTGGTTTGTGGAAGCGG
<i>MsG0180004772.01</i>	CAGCACCACCGTTGAATTGG	GGGCCCATGTAAAGGGCTTA
<i>MsG0180004541.01</i>	CCGCTGAGGTTCTTCCTGTT	ACCAACCAATCGACGACACA
<i>MsG0880047697.01</i>	TTCTTTTTGCGCCACCTTCC	TCAGCCAGGTTCTCAAGTTCC
<i>Ms-ACTIN2</i>	CAAAAGATGGCAGATGCTGAGGAT	CATGACACCAGTATGACGAGGTCG