



Brief Report Evaluation of the Analytical Performance of a Lateral Flow Assay for the Detection of Anti-Coccidioides Antibodies in Human Sera—Argentina

Mariana N. Viale ^{1,*}, Diego H. Caceres ^{2,3,4}, Patricia E. Mansilla ⁵, María C. Lopez-Joffre ¹, Flavia G. Vivot ¹, Andrea N. Motter ⁶, Adriana I. Toranzo ¹ and Cristina E. Canteros ¹

- ¹ Laboratorio Nacional de Referencia en Micología Clínica, Departamento Micología, Instituto Nacional de Enfermedades Infecciosas (INEI), Administración Nacional de Laboratorios e Institutos de Salud (ANLIS) "Dr. Carlos G. Malbrán", Buenos Aires 1282, Argentina; mcecilj@hotmail.com (M.C.L.-J.); fgvivot@anlis.gob.ar (F.G.V.); atoranzo@hotmail.com (A.I.T.); cecanteros@gmail.com (C.E.C.)
- ² Center of Expertise in Mycology Radboudumc/CWZ, 6525GA Nijmegen, The Netherlands; diegocaceres84@gmail.com
- ³ Studies in Translational Microbiology and Emerging Diseases (MICROS) Research Group, School of Medicine and Health Sciences, Universidad del Rosario, Bogota 111221, Colombia
- ⁴ Immuno-Mycologics Inc. (IMMY), Norman, OK 73069, USA
- ⁵ Hospital Interzonal San Juan Bautista, San Fernando del Valle de Catamarca 4700, Argentina; p.evangelina.mansilla@gmail.com
- ⁶ Unidad Operativa Centro de Contención Biológica, Administración Nacional de Laboratorios e Institutos de Salud (ANLIS) "Dr. Carlos G. Malbrán", Buenos Aires 1282, Argentina; andreamotter@gmail.com
- * Correspondence: viale.mariana.noelia@gmail.com; Tel.: +54-01137518090

Abstract: Coccidiomycosis is a potentially life-threatening fungal infection endemic to certain regions of Argentina. The infection is caused by *Coccidioides* spp. and is primarily diagnosed by *Coccidioides* antibody (Ab) detection. Access to rapid, highly accurate diagnostic testing is critical to ensure prompt antifungal therapy. The sona *Coccidioides* Ab Lateral Flow Assay (LFA) performs faster and requires less laboratory infrastructure and equipment compared with other Ab detection assays, potentially providing a substantial improvement for rapid case screening in coccidioidomycosis-endemic regions; however, validation of this test is needed. Thus, we aimed to evaluate the analytical performance of the sona *Coccidioides* Ab (LFA) and compare agreement with anti-*Coccidioides* Ab detection assays. A total of 103 human sera specimens were tested, including 25 specimens from patients with coccidioidomycosis and 78 from patients without coccidioidomycosis. The sona *Coccidioides* Ab Lateral Flow Assay (LFA) was performed with a sensitivity of 88%, and specificity and accuracy of 87%. Furthermore, the *Coccidioides* Ab LFA had good agreement with other anti-*Coccidioides* Ab detection assays. Our findings suggest the sona *Coccidioides* Ab LFA has satisfactory performance and may be useful for diagnosing coccidioidomycosis in endemic regions.

Keywords: coccidioidomycosis; *Coccidioides*; human; antibody; immunodiagnostic; lateral flow assay; enzyme immunoassay; immunodiffusion; counterimmunoelectrophoresis; Argentina

1. Introduction

Coccidioidomycosis is a fungal disease caused by *Coccidioides immitis* and *C. posadasii*. Infection is the result of exposure to contaminated soil, primarily by inhalation of fungal arthroconidia found in the matter. The mycosis predominantly affects individuals who reside in endemic regions, which are localized in the southwestern United States, Northern Mexico, Guatemala, Honduras, Venezuela, Northeastern Brazil, Southern Bolivia, Western Paraguay, and the arid pre-Andes region of Argentina [1–4].

Endemic regions in Argentina are characterized by arid to semi-arid climates, in the arid foothills from the province of Jujuy to the province of Río Negro, with the highest



Citation: Viale, M.N.; Caceres, D.H.; Mansilla, P.E.; Lopez-Joffre, M.C.; Vivot, F.G.; Motter, A.N.; Toranzo, A.I.; Canteros, C.E. Evaluation of the Analytical Performance of a Lateral Flow Assay for the Detection of Anti-*Coccidioides* Antibodies in Human Sera—Argentina. *J. Fungi* **2024**, *10*, 322. https://doi.org/ 10.3390/jof10050322

Academic Editor: Theo N. Kirkland

Received: 19 March 2024 Revised: 17 April 2024 Accepted: 23 April 2024 Published: 28 April 2024



Copyright: © 2024 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https:// creativecommons.org/licenses/by/ 4.0/). number of cases reported in the province of Catamarca [5]. Occupational exposure, particularly during agricultural activities and construction, is associated with a higher risk of infection [3,6–10]. Other activities such as digging, excavation, and disturbance of soil can lead to the inhalation of arthroconidia. Furthermore, climatic variations and environmental changes may influence the distribution of *Coccidioides* [11–14]. Surveillance efforts in Argentina, like the mandatory reporting of cases in endemic regions like Catamarca and the recent incorporation of national-level reporting, aim to monitor epidemiological trends of the disease, identify high-risk areas, and enhance public awareness [15,16].

Reference diagnosis of coccidioidomycosis is performed by microscopic observation of the fungus and culture. Handling of culture specimens is complex due to the high biohazard risk; microscopic analysis is rapid, but the performance of this test is variable. Antibody (Ab) detection using complement fixation (CF), immunodiffusion (ID), counterimmunoelectrophoresis (CIE), and enzyme immunoassay (EIA) is widely adopted, particularly due to the use of minimally invasive specimens and their commendable analytical performance. However, a notable limitation of these assays is the extended turnaround time required to obtain antibody detection results [1,17]. For these reasons, accurate assays are needed. The principal goal of this study was to evaluate the analytical performance of a novel lateral flow assay (LFA) for the rapid detection of anti-*Coccidioides* antibodies in human sera.

2. Materials and Methods

Specimens: 103 human serum specimens were tested. Coccidioidomycosis specimens include 25 sera, these sera came from 12 patients with proven coccidioidomycosis by culture (n = 17 sera), and 8 patients with positive antibody testing for coccidioidomycosis (n = 8 sera). We also tested 78 sera from persons without coccidioidomycosis, including 30 sera specimens from people with other etiologies for respiratory disease who live in endemic regions for coccidioidomycosis. Also, we included 18 sera from patients with positive antibody testing for histoplasmosis, 16 from patients with positive antibody testing for aspergillosis (Figure 1).




Figure 1. Flow chart of specimens tested for the evaluation of the analytical performance of the sona *Coccidioides* Ab Lateral Flow Assay in human sera.

Specimens were remnants from serological studies conducted at the National Reference Laboratory in Clinical Mycology of Argentina (LNRM). These specimens were stored at -20 °C, in single-use aliquots, and were collected between 2009 and 2021. We tested all available sera from patients with proven and positive antibody testing for coccidioidomycosis, and sera from patients with positive antibody testing against other fungi, *Histoplasma, Paracoccidioides*, and *Aspergillus*, as specificity controls. Prior to EIA and LFA testing, specimens were re-tested by ID and CIE methods to validate specimen stability and the presence of anti-*Coccidioides* antibodies. (Figure 1).

Coccidioides Antibody (Ab)-detection Assays: four assays were evaluated. Two inhouse assays, including CIE and ID. These assays were performed using an in-house *Coccidioides* antigen (Ag), and a rabbit serum containing polyclonal Ab for positive control, both produced by the LNRM [18]. Two commercial immunodiagnostic assays from IMMY (Norman, OK, USA) were also used: the Clarus *Coccidioides* Ab Enzyme Immunoassay (*Coccidioides* Ab EIA) and the sona *Coccidioides* Ab Lateral Flow Assay (*Coccidioides* Ab LFA).

Coccidioides Ab CIE was performed following laboratory standard operating procedure (SOP). First, gel plates for CIE were prepared using 1% agarose tipe II medium EEO (Sigma, St Louis, USA)/Tris Acetate-EDTA (TAE) buffer. Wells were created using a puncher and were filled with the specimens (8 μ L), the Ag (2.5 μ L), and the positive control (8 μ L). Electrophoresis was performed in a TAE buffer at 100 V for 90 min. Results were interpreted by the visualization of precipitin bands in the plates.

Coccidioides Ab ID was performed following laboratory SOP. In brief, 20 µL of human sera, *Coccidioides* Ag, and positive control sera were placed into 1% phenolyzed agar ID plates. Plates were incubated at room temperature, and results were interpreted after 72 h. A positive result was interpreted as the presence of precipitin identity bands in the ID plates. Standard operating procedure (SOP) available at http://www.anlis.gov.ar/inei/micologia/?page_id=1000 (accessed on 18 March 2024).

Coccidioides Ab EIA was achieved following the instructions described in the kit package insert. Briefly, 100 µL of diluted 1:441 sera, CF and TP calibrator cutoff, and kit positive and negative controls were placed into microwells coated with CF or TP Coccidioides antigen, respectively. Plates were incubated at 20–25 °C (room temperature) for 30 min. After incubation, plates were manually washed three times using $1 \times$ wash buffer. Then, 100 µL of CF or TP enzyme conjugate (HRP) was added to their respective plates. Plates were incubated for 30 min at room temperature. Following the incubation, plates were washed again three times as previously described. Then, 100 μ L of TMB Substrate was added to each well and incubated in darkness for 10 min at 20–25 °C. Finally, 100 µL of stop solution was dispensed to all wells. Plate's optical densities (OD) were obtained using a dual-wavelength plate reader at 450 and 630 nm using the Labsystems Multiskan RC spectrophotometer (Thermo Electron Corporation, Vantaa, Finland). EIA units were obtained by dividing the OD value of serum and control wells by the blanked OD from the calibrator cutoff well. Specimens with EIA unit values \geq 1.5 were considered positive, specimens with EIA unit values between \geq 1.0 to 1.49 were considered indeterminate, and specimens with EIA unit values <1.0 were considered negative. For the 2 \times 2 table analysis, we considered indeterminate as negatives.

Coccidioides Ab LFA was carried out following manufacturer instructions. Sera specimens were diluted 1:441 using kit specimen diluent. LFA strip was placed in the tube containing 100 μ L of the diluted specimen, and results were interpreted after 30 min of incubation at room temperature. The presence of two pink or red lines, at the control and test zone, was considered a positive result. The presence of only the control line was considered a negative result. The absence of the control line invalidated the result, regarding the presence or the absence of the test line.

Statistical analysis: Test performance and respective 95% confidence intervals (95% CI) were calculated using 2 × 2 tables. Concordance analysis was conducted by calculation of Kappa. Kappa values were interpreted as 0.0 to 0.2, no agreement; 0.21 to 0.39, minimal agreement; 0.40 to 0.59, weak agreement; 0.60 to 0.79, moderate agreement; 0.80 to 0.89, strong agreement; and \geq 0.90, almost perfect agreement. Data were analyzed using STATA 11.

3. Results

Characteristics of patients with coccidioidomycosis: A total of 20 patients with coccidioidomycosis were included, comprising 11 females and 9 males, with ages ranging from 10 to 73 years old. The participants were primarily from Catamarca (n = 12). Regarding baseline conditions, nine patients had no reported conditions, three had HIV, two had diabetes, two had cancer, and four had no available information. Twelve patients presented pulmonary coccidioidomycosis, subacute and chronic forms, five had disseminated coccidioidomycosis, and three had no specified information (Table 1).

Male	<i>n</i> = 9	
Female	n = 11	
Age range	10 to 73 years old	
Residence		
Catamarca	<i>n</i> = 12	
Salta	<i>n</i> = 3	
Santiago del Estero	<i>n</i> = 2	
Mendoza	n = 1	
La Rioja	n = 1	
Jujuy	n = 1	
Baseline conditions		
No reported	<i>n</i> = 9	
HIV	<i>n</i> = 3	
Diabetes	<i>n</i> = 2	
Cancer	<i>n</i> = 2	
No information	n = 4	
Coccidioidomycosis clinical form		
Pulmonary (subacute or chronic)	12 *	
Disseminated	5	
No information	3	

Table 1. Characteristics of patients with coccidioidomycosis.

(*) no acute pulmonary coccidioidomycosis cases recorded.

Coccidioides Ab LFA: Of the 25 sera from patients with coccidioidomycosis, 22 were positive for *Coccidioides* Ab LFA, showing a sensitivity of 88% (95% CI: 69–97). Of the three false negative specimens, one was from a patient with advanced HIV disease; this patient had proven coccidioidomycosis by culture (*C. posadasii*) and also had positive IgM EIA (3.1 EIA units). The second specimen was from a probable case, with positive CIE and IgM EIA (2.4 EIA units). The third specimen was from a culture-proven case (*C. posadasii*), with positive CIE, ID, and IgG EIA (4.1 EIA units). In the 78 controls without coccidioidomycosis, 10 specimens tested positive by the LFA, showing a specificity of 87% (95% CI: 78–94). Most of the false positive results were found in sera from people with paracoccidioidomycosis, n = 7 (7 out of 16 sera, 44% cross-reactivity), one false positive in a patient with histoplasmosis (1 out of 18 sera, 6% cross-reactivity), and one false positive in sera from a person with respiratory disease who live in a endemic region for coccidioidomycosis (1 out of 30, 3% false positivity). Assay accuracy was 87% (95% CI: 79–93) (Table 2A, and Figure 2).

We also evaluate the performance of the other assays for the anti-*Coccidioides* Ab detection. The performance of these assays is described below.

Coccidioides Ab CIE: In total, 23 of the 25 sera from individuals with coccidioidomycosis had positive results for CIE (sensitivity 92%, 95% CI: 74–99). The two false negatives were proven cases. One serum was a proven case that tested positive by the LFA only. The second specimen was a proven case with positive IgM EIA only (HIV patient described above). On the other hand, 77 out of the 78 non-coccidioidomycosis control specimens were negative (99% specificity, 95% CI: 93–100). The only false positive result was found in a patient with proven paracoccidioidomycosis. The CIE shows an accuracy of 97% (95% CI: 92–99) (Table 2B, and Figure 2).

A) Coccidioid	es Ab LFA				
			Coccidioidomycosis		
			+	_	
	LFA	+	22	10	
	LFA	_	3	68	
tatistic			Value (%)	95% CI	
ensitivity			88	69–97	
pecificity			87	78–94	
ccuracy			87	79–93	
3) Coccidioid	es Ab CIE				
			Coccidioidomycosis		
			+	_	
	CIE	+	23	1	
	CIE	_	2	77	
atistic			Value (%)	95% CI	
ensitivity			92	74–99	
pecificity			99	93–100	
curacy			97	92–99	
) Coccidioid	es Ab ID				
			Coccidioidomycosis		
			+	_	
	ID	+	17	0	
	ID.	_	8	78	
tistic			Value (%)	95% CI	
sitivity			68	47-85	
cificity			100	95–100	
curacy			92	85–97	
Coccidioid	les Ab EIA (Ig	;M+IgG)			
			Coccidioidomycosis		
			+	_	
	EIA	+	24	2	
		_	1	76	
atistic			Value (%)	95% CI	
ensitivity			96	80–100	
pecificity			97	91–100	
ccuracy			97	92–99	

Table 2. Analytical performance of assays for detection of Coccidioides Ab in human sera.

Coccidioides Ag ID: In total, 17 out of 25 sera tested positive (sensitivity 68%, 95% CI: 47-85). False negatives involved six specimens from proven cases, and two specimens from probable cases. Four sera from proven cases that tested positive for CIE, LFA, and IgG EIA. One serum from a probable case that tested positive by CIE, LFA, IgM, and IgG EIA. One specimen from a proven case tested positive by the IgM EIA only (patient with HIV described before). And one serum from a proven case that tested positive only by the LFA. No false positive results were recorded using ID (specificity of 100%; 95% CI: 95–100). ID displayed an assay accuracy was 92% (95% CI: 85–97) (Table 2C, and Figure 2).






Figure 2. Comparison of turn-around of results and sensitivity, specificity, and accuracy of assay for the detection of *Coccidioides* Ab.

Coccidioides Ab EIA: IgM EIA was performed with 36% sensitivity, 9 positive out of 25 coccidioidomycosis cases. IgG was performed with a sensitivity of 88%, with 22 positives out of 25 cases (the 3 discrepant results were indeterminate). Combining IgM and IgG, performance was 96%, showing one false negative result. This IgM/IgG false negative was serum from a proven case, with indeterminate IgG (1.0 EIA units), negative IgM, positive LFA, and negative CIE and ID. In terms of assay specificity, two false positives were observed using IgM EIA (IgM specificity of 97%; 95% CI: 91–100), these false positives were observed using IgG EIA. Combining IgM/IgG EIA, we found a specificity of 97% (95% CI: 91–100). This specificity was impacted by the histoplasmosis cross-reactions observed using IgM (Table 2D, and Figure 2). It is important to highlight that in control specimens the presence of indeterminate results was frequent. We observed 7 sera with IgM indeterminate result and 12 sera with IgG indeterminate result. However, the ROC analysis did not show an alternative cut-off that improved EIA performance, the reason why we kept the cut-off recommended by the manufacturer.

Concordance analysis: The concordance between all assays ranged from moderate to very good. The majority of assays were performed with good agreement (Figure 3).





Most of the coccidioidomycosis cases are diagnosed using Ab detection assays, based on complement fixation (CF), gel precipitation tests (CIE and ID), and EIA. The main limitation of these assays is the need for highly trained laboratorians, special laboratory instruments, and prolonged turnaround time for results [19]. LFA technology has significantly improved this limitation; these assays are simple, accurate, and rapid. In this study, we were able to evaluate and compare a LFA for the detection of anti-*Coccidioides* Ab in human sera, sōna *Coccidioides* Ab LFA (IMMYTM, Norman, OK, USA). Results here showed high analytical performance and good correlation with other assays used for the immunodiagnosis of coccidioidomycosis.

The sona *Coccidioides* Ab LFA was performed with a sensitivity of 88%. In terms of sensitivity, these results are comparable to results reported by Contreas et al. at the University of California, in Los Angeles, CA, USA (poster ASM meeting 2019), but are substantially different from results reported by Donovan et al. at the University of Arizona, AZ, USA [20,21]. In the present study, it is important to highlight that in two of three false negative LFA specimens, we observed positive IgM results and negative IgG results, suggesting early disease in false negatives. Another factor could be the high levels of immunosuppression, as observed in patients with advanced HIV disease. These factors might have affected the concentration of anti-*Coccidioides* Ab, which was likely below the limit of detection of the LFA [22]. In addition, the main substantial difference of this study compared with the report of Donovan et al. is the low frequency of patients with chronic and disseminated coccidioidomycosis. That could explain the higher LFA sensitivity in this study.

We also evaluated the performance of the sensitivity of other immunodiagnostic assays. Based on that analysis, we found that the combination of IgM and IgG detection was the assay with the highest sensitivity, followed by the CIE. ID was 68% sensitivity, and the detection of IgM alone was the assay with the lowest sensitivity (36%).

Regarding specificity, we observed that cross-reactivity with sera from patients with paracoccidioidomycosis was the most frequent discrepancy. We also observed a lower frequency of false positive results in sera from patients with aspergillosis and histoplasmosis. It is well known that cross-reactions, especially in specimens from individuals with fungal infections, could happen. One false positive result was observed in serum from a patient with non-coccidioidomycosis respiratory disease, this patient resided in a endemic region for coccidioidomycosis, probably the positive results was a consequence of the environmental exposure to *Coccidioides* product and of living in an endemic region. By eye reading, differences in band intensity are difficult to interpret, a reason why the addition of a LFA reader or other alternatives for automatic reading of LFA devices could improve the interpretation and performance of the assay.

Simplicity, rapid testing, and high accuracy are characteristics that make the sōna *Coccidioides* Ab Lateral Flow Assay a desirable alternative for rapid screening of possible coccidioidomycosis cases. But, compared with the other Ab assays evaluated in this study, the *Coccidioides* Ab LFA performed with the lowest specificity. For this reason, it would be important to confirm positive LFA results by the use of a confirmatory assay, such as the *Coccidioides* Ab EIA, or the *Coccidioides* Ab CIE, given the high diagnostic analytical performance of these assays. This finding was previously reported, and similar recommendations have been proposed [23].

Principal limitations of this study were the small sample size, the use of sera from some probable cases, the limited access to charts of cases, the unblinded testing of specimens, the lack of an automatic reader for LFA interpretation, and the lack of sera specimens of asymptomatic people who are residents and not residents in the coccidioidomycosisendemic region. Another major limitation is the lack of specimens from patients with acute pulmonary coccidioidomycosis.

5. Conclusions

Here, we evaluate a novel LFA for the detection of anti-*Coccidioides* antibodies and were able to evaluate the analytical performance of other assays for the immunodiagnosis of coccidioidomycosis. Diagnosis of coccidioidomycosis could take days to hours using CF, gel precipitation tests, and EIA; the use of LFAs reduces this time to minutes. Another advantage of LFAs is the simplicity of testing implementation in low-complexity laboratories located in regions where the disease is endemic, with limited access to reference laboratories.

Author Contributions: Conceptualization, M.N.V., D.H.C. and C.E.C.; Methodology, M.N.V., D.H.C., P.E.M., F.G.V., A.N.M., A.I.T. and C.E.C.; Software, M.N.V., D.H.C., P.E.M. and C.E.C.; Validation, M.N.V., D.H.C., P.E.M., M.C.L.-J., F.G.V., A.N.M., A.I.T. and C.E.C.; Formal Analysis, M.N.V., D.H.C., P.E.M., M.C.L.-J., F.G.V., A.N.M., A.I.T. and C.E.C.; Investigation, M.N.V., D.H.C., P.E.M., M.C.L.-J., F.G.V., A.N.M., A.I.T. and C.E.C.; Presources, M.N.V., D.H.C., P.E.M. and C.E.C.; Data Curation, M.N.V., D.H.C., P.E.M. and C.E.C.; Writing—Original Draft, M.N.V., D.H.C., P.E.M. and C.E.C.; Writing—Review and Editing, M.N.V., D.H.C., P.E.M., M.C.L.-J., F.G.V., A.N.V., D.H.C., P.E.M. and C.E.C.; Visualization, M.N.V., D.H.C., P.E.M. and C.E.C.; Visualization, M.N.V., D.H.C., P.E.M. and C.E.C.; Funding Acquisition, M.N.V., D.H.C. and C.E.C. All authors have read and agreed to the published version of the manuscript.

Funding: This work was partially supported by the project FOCANLIS 2019 NRU 1940 "Vigilancia activa de la Coccidioidomicosis en la principal área endémica de Argentina a través de la optimización del diagnóstico y de las estrategias de búsqueda de fuentes de infección".

Institutional Review Board Statement: The study was conducted in accordance with the Declaration of Helsinki, and approved by the Institutional Review Board of the Instituto de Epidemiología (INE), Administración Nacional de Laboratorios e In-stitutos de Salud "Dr. Carlos G. Malbrán" (ANLIS).

Informed Consent Statement: Specimens were remnants from serological studies conducted at the National Reference Laboratory in Clinical Mycology of Argentina. All samples were coded and anonymized, so they were unlinked from any sensitive data from patients.

Data Availability Statement: Data are contained within the article.

Acknowledgments: We want to thank the members of the National Mycology Laboratory Network of Argentina. We also want to thank Jeremy A. W. Gold for the critical review and editing of the manuscript.

Conflicts of Interest: Author Diego H. Caceres was employed by the company IMMY since November 2021 (Latin America account manager). The remaining authors declare that the research was conducted in the absence of any commercial or financial relationships that could be construed as a potential conflict of interest.

References

- 1. Williams, S.L.; Chiller, T. Update on the Epidemiology, Diagnosis, and Treatment of Coccidioidomycosis. *J. Fungi* 2022, *8*, 666. [CrossRef] [PubMed]
- 2. Bays, D.J.; Thompson, G.R., 3rd. Coccidioidomycosis. Infect. Dis. Clin. N. Am. 2021, 35, 453–469. [CrossRef] [PubMed]
- Lockhart, S.R.; Toda, M.; Benedict, K.; Caceres, D.H.; Litvintseva, A.P. Endemic and Other Dimorphic Mycoses in The Americas. J. Fungi 2021, 7, 151. [CrossRef] [PubMed]
- 4. Thompson, G.R.; Pasqualotto, A.C. Endemic mycoses: Expansion of traditional geographic ranges and pitfalls in management. *Mycoses* **2021**, *64*, 989–992. [CrossRef] [PubMed]
- 5. Canteros, C.; Toranzo, A.; Ibarra-Camou, B.; David, V.; Carrizo, S.G.; Santillán-Iturres, A.; Serrano, J.; Fernández, N.; Capece, P.; Gorostiaga, J.; et al. Coccidioidomycosis in Argentina, 1892–2009. *Rev. Argent. Microbiol.* **2010**, *42*, 261–268. [PubMed]
- 6. Smith, D.J.; Free, R.J.; Thompson, G.R., 3rd; Baddley, J.W.; Pappas, P.G.; Benedict, K.; Gold, J.A.W.; Andes, D.; Bradsher, R.W.; Ebell, M.H.; et al. Clinical Testing Guidance for Coccidioidomycosis, Histoplasmosis, and Blastomycosis in Patients with Community-Acquired Pneumonia for Primary and Urgent Care Providers. *Clin. Infect. Dis.* 2023. *ahead of print*. [CrossRef]
- de Perio, M.A.; Benedict, K.; Williams, S.L.; Niemeier-Walsh, C.; Green, B.J.; Coffey, C.; Di Giuseppe, M.; Toda, M.; Park, J.-H.; Bailey, R.L.; et al. Occupational Histoplasmosis: Epidemiology and Prevention Measures. *J. Fungi* 2021, 7, 510. [CrossRef] [PubMed]

- Toda, M.; Gonzalez, F.J.; Fonseca-Ford, M.; Franklin, P.; Huntington-Frazier, M.; Gutelius, B.; Kawakami, V.; Lunquest, K.; McCracken, S.; Moser, K.; et al. Notes from the Field: Multistate Coccidioidomycosis Outbreak in U.S. Residents Returning from Community Service Trips to Baja California, Mexico—July–August 2018. MMWR Morb. Mortal. Wkly. Rep. 2019, 68, 332–333. [CrossRef] [PubMed]
- Caceres, D.H.; Tap, R.M.; Alastruey-Izquierdo, A.; Hagen, F. Detection and Control of Fungal Outbreaks. *Mycopathologia* 2020, 185, 741–745. [CrossRef] [PubMed]
- Viale, M.N.; López-Joffre, M.C.; Motter, A.N.; Mansilla, P.E.; Vivot, F.G.; Acevedo, F.M.; del Valle David, V.; Carrizo, S.G.; Serrano, J.; Chacón, Y.; et al. Epidemiology of coccidioidomycosis in Argentina, an update. *Med. Mycol.* 2024, 2024, myae024. [CrossRef] [PubMed]
- 11. Comrie, A.C. Climate factors influencing coccidioidomycosis seasonality and outbreaks. *Environ. Health Perspect.* **2005**, *113*, 688–692. [CrossRef] [PubMed]
- Brown, J.; Benedict, K.; Park, B.J.; Thompson, G.R., 3rd. Coccidioidomycosis: Epidemiology. *Clin. Epidemiol.* 2013, 5, 185–197. [PubMed]
- 13. Gorris, M.E.; Neumann, J.E.; Kinney, P.L.; Sheahan, M.; Sarofim, M.C. Economic Valuation of Coccidioidomycosis (Valley Fever) Projections in the United States in Response to Climate Change. *Weather Clim. Soc.* **2021**, *13*, 107–123. [CrossRef] [PubMed]
- Park, B.J.; Sigel, K.; Vaz, V.; Komatsu, K.; McRill, C.; Phelan, M.; Colman, T.; Comrie, A.C.; Warnock, D.W.; Galgiani, J.N.; et al. An epidemic of coccidioidomycosis in arizona associated with climatic changes, 1998–2001. J. Infect. Dis. 2005, 191, 1981–1987. [CrossRef] [PubMed]
- de SA, Información Jurídica: Programa Provincial de Prevención DyTdIC. LEY 5.523. SAN FERNANDO DEL VALLE DE CATAMARCA, 23 de Noviembre de 2017. Boletín Oficial, 9 de Enero de 2018. Vigente, de Alcance General. Id SAIJ: LPK0005523 2017. Available online: http://www.saij.gob.ar/5523-local-catamarca-programa-provincial-prevencion-deteccion-tratamientococcidioidomicosis-lpk0005523-2017-11-23/123456789-0abc-defg-325-5000kvorpyel?utm_source=newsletter-mensual&utm_ medium=email&utm_term=mensual&utm_campaign=ley-provincial (accessed on 18 March 2024).
- 16. Ministerio de Salud RA. Manual de Normas Yprocedimientos de Vigilancia y Control de Eventos de Notificación Obligatoria, Ministerio de Salud Argentina 2022. Available online: https://bancos.salud.gob.ar/recurso/manual-de-normas-y-procedimientosde-vigilancia-y-control-de-eventos-de-notificacion (accessed on 18 March 2024).
- 17. Saubolle, M.A. Laboratory Aspects in the Diagnosis of Coccidioidomycosis. *Ann. N. Y. Acad. Sci.* **2007**, 1111, 301–314. [CrossRef] [PubMed]
- 18. Perrotta, D.; Vivot, W.; Lee, W.; Rivas, C.; Yabo, M.; Rodero, L.; Canteros, C.; Davel, G. Specific fungal antisera production in rabbits. *Rev. Argent. Microbiol.* **1998**, *30*, 115–121. [PubMed]
- McHardy, I.H.; Barker, B.; Thompson, G.R., 3rd. Review of Clinical and Laboratory Diagnostics for Coccidioidomycosis. J. Clin. Microbiol. 2023, 61, e0158122. [CrossRef] [PubMed]
- Contreras, D.A.; Li, Q.; Garner, O.B. Evaluation of Sona Lateral Flow Assay for the Rapid Detection of Coccidioides Immitis; ASM Microbe: San Francisco, CA, USA, 2019.
- Donovan, F.M.; Ramadan, F.; Khan, S.; Bhaskara, A.; Lainhart, W.D.; Narang, A.T.; Mosier, J.M.; Ellingson, K.D.; Bedrick, E.J.; Saubolle, M.; et al. Comparison of a Novel Rapid Lateral Flow Assay to Enzyme Immunoassay Results for Early Diagnosis of Coccidioidomycosis. *Clin. Infect. Dis.* 2020, 73, e2746–e2753. [CrossRef] [PubMed]
- 22. Caceres, D.H.; Chiller, T.; Lindsley, M.D. Immunodiagnostic Assays for the Investigation of Fungal Outbreaks. *Mycopathologia* **2020**, *185*, 867–880. [CrossRef] [PubMed]
- 23. Lu, J.; Banaei, N. Investigation of discordant positive Sona lateral flow assay results for detection of *Coccidioides* antibodies. *J. Clin. Microbiol.* **2023**, *61*, e0062323. [CrossRef]

Disclaimer/Publisher's Note: The statements, opinions and data contained in all publications are solely those of the individual author(s) and contributor(s) and not of MDPI and/or the editor(s). MDPI and/or the editor(s) disclaim responsibility for any injury to people or property resulting from any ideas, methods, instructions or products referred to in the content.