

# Synthesis and improved cross-linking properties of C5-modified furan bearing PNAs

Joke Elskens<sup>1</sup>, Alex Manicardi<sup>1,2</sup>, Valentina Costi<sup>2</sup>, Annemieke Madder<sup>1\*</sup>, Roberto Corradini<sup>2\*</sup>

<sup>1</sup> Organic and Biomimetic Chemistry Research Group, Department of Organic and Macromolecular Chemistry, Ghent University, Krijgslaan 281-S4, Gent 9000, Belgium.

<sup>2</sup> Department of Chemistry, Life Sciences and Environmental Sustainability, University of Parma, Parco Area delle Scienze 17/A, Parma 43124, Italy

## Supplementary Material - Table of contents

<b>Solid Phase Diels-Alder optimization</b>	<b>2-5</b>
Solvent and maleimide concentration optimization	2
Figure S1: UPLC analysis of products obtained using N-(N-Boc-2-aminoethyl)maleimide as protecting group	3
Maleimide concentration and equivalents optimization	3
Figure S2: relative by-product formation	4
Maleimide equivalents and reaction time optimization	4
Figure S3: relative by-product formation	5
Visualization of the retro-DA reaction time course	6
Figure S4: time course of the retro-DA reaction	6
<b>HPLC-DAD-HRMS chromatograms</b>	<b>7</b>
Figure S5: Overview of the PNA sequences and furan-PNA building blocks	7
Figure S6: HPLC-DAD-HRMS chromatogram of purified <b>PNA 2</b>	8
Figure S7: HPLC-DAD-HRMS chromatogram of purified <b>PNA 3</b>	9
Figure S8: HPLC-DAD-HRMS chromatogram of purified <b>PNA 5</b>	10
Figure S9: HPLC-DAD-HRMS chromatogram of purified <b>PNA 6</b>	11
Figure S10: HPLC-DAD-HRMS chromatogram of purified <b>PNA 8</b>	12
Figure S11: HPLC-DAD-HRMS chromatogram of purified <b>PNA 9</b>	13
<b>PAGE analysis</b>	<b>14</b>
Figure S12: PAGE analysis of the cross-link experiment between <b>ON1-ON4</b> and <b>PNA1-PNA3</b>	14
Figure S13: PAGE analysis of the cross-link experiment between <b>ON1-ON4</b> and <b>PNA4-PNA6</b>	15
Figure S14: PAGE analysis of the cross-link experiment between <b>ON1-ON4</b> and <b>PNA7-PNA9</b>	15
Figure S15: PAGE analysis of the strand displacement experiment between <b>ON1-ON2</b> and <b>PNA1-PNA9</b>	16
<b>MALDI analysis</b>	<b>17</b>
Figure S16: MALDI analysis of the cross-link product of <b>ON2</b> and <b>PNA4</b>	17
Figure S17: HPLC trace and MALDI analysis of the cross-link product of <b>ON1</b> and <b>PNA2</b>	18
<b>Mechanism of Interstrand Crosslink formation</b>	<b>19</b>
Figure S18: Schematic representation of interstrand crosslink formation	19

### Solid Phase Diels-Alder optimization

In order to optimize the Diels-Alder reaction conditions on the solid phase, a small PNA sequence was N-capped with 2-furanpropionic acid (Fur-GATCT-Gly-NH<sub>2</sub>, loading 7,45 mg/μmol).

### Solvent and maleimide concentration optimization

#### General protocol

- 1- Load a 200 μL Eppendorf with 0.2 μmol resin (approx. 1.86 mg);

- 2- Prepare desired solutions of N-(N-Boc-2-aminoethyl)maleimide, and transfer them to the Eppendorf tube;
- 3- Allow Diels-Alder reaction for 15 hours at 90°C
- 4- Transfer the resin beads in a 1 mL SPE tube equipped with PE frit and wash the resin with DMF/DCM;
- 5- Add 100 µL of a 10% m-cresol solution in TFA and allow to react for 1h30';
- 6- Percolate the solution in a 2 mL Eppendorf tube, and wash the resin with 80 µL TFA;
- 7- Add 1.8 mL ethyl ether and allow to precipitate for 2h at -20 °C;
- 8- Recover the precipitate by centrifugation, and wash the pellet with ethyl ether;
- 9- Redissolve the dry pellet in 500 µL mQ, transfer in a 500 µL Eppendorf tube and heat for 15 hours at 90 °C;
- 10- Submit the samples to UPLC-MS analysis.


#### Solution used

Eq. maleimide / Concentration	Solvent: Toluene	Solvent: DMF
100 eq / 2 M	Test A1	Test B1
50 eq / 1 M	Test A2	Test B2
20 eq / 0.5 M	Test A3	Test B3
10 eq / 0.25 M	Test A4	Test B4
5 eq / 0.125 M	\	Test B5

#### Target identification

PNA \ Charge	No charge	2+	3+	4+
Target	1546.5	774	516	387
Alkylated	1712.7	857	571	429
Cycloaddition by-prod	1686.6	844	563	422

Results: when toluene is used as a solvent, it is not possible to prevent furan alkylation or the formation of a by-product with an additional unit of maleimide connected to the PNA strand. When using DMF it is possible to prevent the alkylation (20 eq / 0.5 M seems enough) but not the formation of the extra by-product. In all cases the alkylation of the furan is inversely proportional to the amount of maleimide used, while the formation of the by-product is directly proportional.



**Figure S1.** UPLC analysis of products obtained for Fur-GATCT-Gly-NH<sub>2</sub> PNA synthesis obtained under various conditions (A1–A4 and B1–B5 tests reported above) using *N*-(*N*-Boc-2-aminoethyl)maleimide as protecting group. Left column: tests A, reference; Right column: tests B. Top panels: UPLC-MS traces; middle panels: MS spectra of 2.40–3.30 minutes region; bottom panels: MS spectra of 5.00–6.00 minutes region. MW target (red boxes) = 1546.5; [MH<sub>2</sub>]<sup>2+</sup> = 774.3; [MH<sub>3</sub>]<sup>3+</sup> = 516.5; MW alkylated product = 1712.7 (blue boxes); [MH<sub>2</sub>]<sup>2+</sup> = 857.4; [MH<sub>3</sub>]<sup>3+</sup> = 571.9 MW cycloadduct = 1686.6 (green boxes); [MH<sub>2</sub>]<sup>2+</sup> = 844.3; [MH<sub>3</sub>]<sup>3+</sup> = 563.2.

## Maleimide concentration and equivalents optimization

### General protocol

- 1- Load a 200  $\mu$ L Eppendorf with 0.2  $\mu$ mol resin (approx. 1.86 mg);
- 2- Prepare desired solutions of *N*-(*N*-Boc-2-aminoethyl)maleimide, and transfer them to the Eppendorf tube;
- 3- Run each condition in triplicate;

- 4- Allow Diels-Alder reaction to proceed for 7.5 hours at 90°C
- 5- Transfer the resin beads in a 1 mL SPE tube equipped with PE frit and wash the resin with DMF/DCM;
- 6- Add 100  $\mu$ L of a 10% m-cresol solution in TFA and allow to react for 1h 30 min;
- 7- Percolate the solution in a 2 mL Eppendorf tube, and wash the resin with 80  $\mu$ L TFA;
- 8- Add 1.8 mL ethyl ether and allow to precipitate for 2h at -20 °C;
- 9- Recover the precipitate by centrifugation, and wash the pellet with ethyl ether;
- 10- Redissolve the dry pellet in 500  $\mu$ L mQ, transfer in a 500  $\mu$ L Eppendorf tube and heat for 5 hours at 90 °C;
- 11- Submit the samples to UPLC-MS analysis.
- 12- Data analysis


#### Solution used

Maleimide concentration	30 eq	20 eq	10 eq
0.5 M	\	Test B1	Test C1
0.25 M	Test A2	Test B2	Test C2
0.1 M	Test A3	Test B3	\

#### Target identification

PNA \ Charge	No charge	2+	3+	4+
Target	1546.5	774	516	387
Alkylated	1712.7	857	571	429
Cycloaddition by-prod	1686.6	844	563	422

**Data analysis:** For each single UPLC-MS trace three different XIC were obtained extracting the signal relative to the 2+ and 3+ m/z signals. Peaks were then integrated using MassLynx 4.0 algorithm. **Results:** 20 equivalent of maleimide at a concentration of 0.25 M gives the best results in term of by-product formation.



**Figure S2.** Relative by-product formation for Fur-GATCT-Gly-NH<sub>2</sub> PNA synthesis obtained under various conditions (A1–A4 and B1–B5 tests reported above) using *N*-(*N*-Boc-2-aminoethyl)maleimide as protecting group. Quantities are normalized to target product intensity in the UPLC-MS analysis (carried out as in figure S1). Data presented are an average of 3 different experiments.

## Maleimide equivalents and reaction time optimization

### General protocol

- 1- Load a 200  $\mu$ L Eppendorf with 0.2  $\mu$ mol resin (approx. 1.86 mg);
- 2- Prepare the desired 0.25 M solution of N-(N-Boc-2-amminoethyl)maleimide, and transfer the correct amount to the Eppendorf tube;
- 3- Run each condition in triplicate;
- 4- Allow Diels-Alder reaction to proceed for the specified time at 90°C
- 5- Transfer the resin beads in a 1 mL SPE tube equipped with PE frit and wash the resin with DMF/DCM;
- 6- Add 100  $\mu$ L of a 10% m-cresol solution in TFA and allow to react for 1h30';
- 7- Percolate the solution in a 2 mL Eppendorf tube, and wash the resin with 80  $\mu$ L TFA;
- 8- Add 1.8 mL ethyl ether and allow to precipitate for 2h at -20 °C;
- 9- Recover the precipitate by centrifugation, and wash the pellet with ethyl ether;
- 10- Redissolve the dry pellet in 500  $\mu$ L mQ, transfer in a 500  $\mu$ L Eppendorf tube and heat for 5 hours at 90 °C;
- 11- Submit the samples to UPLC-MS analysis.
- 12- Data analysis

### Conditions used


Maleimide amount	15 h at 90°C	5 h at 90 °C
10 eq	Test A1	Test B1
20 eq	Test A2	Test B2
30 eq	Test A3	\

### Target identification

PNA \ Charge	No charge	2+	3+	4+
Target	1546.5	774	516	387
Alkylated	1712.7	857	571	429
Cycloaddition by-prod	1686.6	844	563	422


**Data analysis:** For each single UPLC-MS trace three different XIC were obtained extracting the signal relative to the 2+ and 3+ m/z signals. Peaks were then integrated using MassLynx 4.0 algorithm.

**Results:** The use of 10 equivalents of maleimide at a concentration of 0.25 M gives the best results in terms of by-product formation, if the Diels-Alder reaction is conducted for a reduced time (5h).



**Figure S3.** Relative by-product formation for Fur-GATCT-Gly-NH<sub>2</sub> PNA synthesis obtained under various conditions (A1–A4 and B1–B5 tests reported above) using *N*-(*N*-Boc-2-aminoethyl)maleimide as protecting group. Quantities are normalized to target product intensity in the UPLC-MS analysis (carried out as in figure S1). Data presented are an average of 3 different experiments.

#### Visualization of retro-DA reaction time course



**Figure S4.** Time-course of the retro-DA reaction. UPLC-MS chromatograms (XIC with *m/z* selected for the Diels alder product (2.77 min) and for the Retro-Diels-Alder product (3.28 min) for reactions carried out for Fur-GATCT-Gly-NH<sub>2</sub> PNA synthesis obtained using 20 equivalents of *N*-(*N*-Boc-2-aminoethyl)maleimide in the protection step (7.5 h at 90°C), and different times of the retro-DA reaction at 90°C after cleavage. a) before retro-DA; b-e) after 1, 3, 4, and 5 h respectively.

## HPLC-DAD-HRMS chromatograms







Figure S5. Overview of PNA sequences and furan-PNA monomers.

RT: 0.00 - 50.00


VC1-75 B #775-978 RT: 13.38-16.30 AV: 102 NL: 1.59E6  
T: FTMS + p ESI Full ms [300.00-2000.00]

**Figure S6.** HPLC-DAD-HRMS chromatogram of purified PNA 2. HPLC-UV at 254 nm trace (top), MS spectrum of the corresponding peak at 15.58 min (center) and mathematical deconvolution of the multicharged signals (insert). (MW: 3158.2976; Molecular formula:  $C_{132}H_{169}N_{69}O_{36}$ ).

RT: 0.00 - 50.00





VC1-50-1B #890-979 RT: 15.40-16.91 AV: 45 NL: 1.45E5  
F: FTMS + p ESI Full ms [300.00-2000.00]



**Figure S7.** HPLC-DAD-HRMS chromatogram of purified PNA 3. HPLC-UV at 254 nm trace (top), MS spectrum of the corresponding peak at 15.93 min (center) and mathematical deconvolution of the multicharged signals (insert). (MW: 3186.3051; Molecular formula:  $C_{132}H_{169}N_{71}O_{36}$ ).





RT: 0.00 - 50.00

VC1-32 #831-1022 RT: 14.36-17.10 AV: 96 NL: 2.68E6  
T: FTMS + p ESI Full ms [300.00-2000.00]



**Figure S9.** HPLC-DAD-HRMS chromatogram of purified PNA 6. HPLC-UV at 254 nm trace (top), MS spectrum of the corresponding peak at 16.03 min (center) and mathematical deconvolution of the multicharged signals (insert). (MW: 3324.3570; Molecular formula: C<sub>127</sub>H<sub>163</sub>N<sub>67</sub>O<sub>35</sub>).

RT: 0.00 - 50.00

VC1-77 #801-1027 RT: 13.83-16.67 AV: 114 NL: 5.65E6  
T: FTMS + p ESI Full ms [300.00-2000.00]

**Figure S10.** HPLC-DAD-HRMS chromatogram of purified PNA 8. HPLC-UV at 254 nm trace (top), MS spectrum of the corresponding peak at 15.60 min (center) and mathematical deconvolution of the multicharged signals (insert). (MW: 3420.3795; Molecular formula: C<sub>137</sub>H<sub>173</sub>N<sub>71</sub>O<sub>38</sub>).

RT: 0.00 - 50.00


VC2-01 #785-1028 RT: 13.56-16.83 AV: 122 NL: 3.83E6  
T: FTMS + p ESI Full ms [300.00-2000.00]

**Figure S11.** HPLC-DAD-HRMS chromatogram of purified PNA 9. HPLC-UV at 254 nm trace (top), MS spectrum of the corresponding peak at 15.77 min (center) and mathematical deconvolution of the multicharged signals (insert). (MW: 3448.3870; Molecular formula:  $C_{137}H_{173}N_{73}O_{38}$ ).


### PAGE analysis

Crosslink samples were analyzed on a 20% polyacrylamide gel (acrylamide:bisacrylamide 19:1) prepared in 1x Tris-Borat-EDTA (TBE) buffer containing 7 M urea. The temperature of the gel was stabilized with a Julabo F12 at 25°C. The power supply used for gel electrophoresis was a consort EV202 and a constant voltage of 230 V was used to run the gels. Gels were stained with SYBR gold (Thermo Fisher Scientific, Life Technologies) and pictures were taken with an Autochemi imaging system (UVP). 4 µL of the crosslink solution (10 µM) were mixed with 16 µL formamide and from this mixture 8 µL were loaded on the gel. It should be noted that staining of PNA is considerably more difficult and much less sensitive than DNA staining. As a result individual ssPNA strands cannot be visualized on PAGE and can thus not be seen in the figures below. Consequently also PNA-DNA crosslinked complexes (boxes in purple) are only visualized with low intensity.


A DNA impurity occurring in all samples appears in all gels as a slower running spot between the cross-linked complex and the ssDNA spots, due to the considerably higher sensitivity of DNA for SybrGold staining




**Figure S12.** PAGE analysis: Cross-link experiment between ON 1-4 and PNA 1-3. Analysis was done on a 20% polyacrylamide gel (acrylamide:bisacrylamide 19:1) prepared in a TBE buffer containing 7M urea. A constant voltage of 230V was applied at a temperature of 25 °C. Gels were stained with a SYBR gold solution and gels were analysed using an Autochemi Imaging system.



**Figure S13.** PAGE analysis: cross-link experiment between ON 1-4 and PNA 4-6. Analysis was done on a 20% polyacrylamide gel (acrylamide:bisacrylamide 19:1) prepared in a TBE buffer containing 7M urea. A constant voltage of 230V was applied at a temperature of 25 °C. Gels were stained with a SYBR gold solution and gels were analysed using an Autochemi Imaging system.




**Figure S14.** PAGE analysis: cross-link experiment between ON 1-4 and PNA 7-9. Analysis was done on a 20% polyacrylamide gel (acrylamide:bisacrylamide 19:1) prepared in a TBE buffer containing 7M urea. A constant voltage of 230V was applied at a temperature of 25 °C. Gels were stained with a SYBR gold solution and gels were analysed using an Autochemi Imaging system.




**Figure S15.** PAGE analysis. **A)** Strand displacement experiment of ON1 and PNA 1-9. **B)** Strand displacement experiment between ON 2 and PNA 1-9. Analysis was done on a 20% polyacrylamide gel (acrylamide:bisacrylamide 19:1) prepared in a TBE buffer containing 7M urea. A constant voltage of 230V was applied at a temperature of 25 °C. Gels were stained with a SYBR gold solution and gels were analysed using an Autochemi Imaging system.

#### MALDI analysis

Prior to MALDI analysis, the crosslinked product was purified using reverse phase HPLC (RP-HPLC). Purification was performed on an Agilent 1200b System equipped with a Waters X-bridge 130 Å Oligonucleotide C18 column (2.5  $\mu$ M, 4.6 mm x 50 mm) at either a column temperature of 50 °C. Acetonitrile and 0.1 M TEAA buffer with 5% acetonitrile were used as mobile phase and were applied through a gradient of 0-18% MeCN in 20 minutes. The isolated crosslinked samples was dissolved in 5  $\mu$ L milliQ. Sample and matrix (2,5-dihydroxybenzoic acid; 15 mg in 100  $\mu$ L milliQ and 50  $\mu$ L MeCN with 1% TFA) were spotted in a 1:1 ratio onto the MALDI plate. The PNA-DNA crosslinked samples were analyzed by MALDI-TOF on an ABI Voyager DE-STR system equipped with a high performance nitrogen laser (337 nm). PNA-DNA crosslinked samples were analyzed in linear, positive mode.




**Figure S16.** MALDI analysis of the cross-link product of **ON 2** and **PNA 4**. The labeled peak is attributed to the DNA-PNA cross-linked duplex (Expected mass: 6591 Da).



**Figure S17.** HPLC trace of the ICL experiment (top panel) and MALDI analysis (bottom panel) of the cross-link product of ON 1 and PNA 2 (peak at 11.00–11.50 min). The labeled peak is attributed to the DNA-PNA cross-linked duplex (Expected mass: 6477 Da). Peak at 8.5 min correspond to unreacted ON 1.

## Mechanism of Interstrand Crosslink formation



**Figure S18.** Schematic representation of interstrand crosslink formation between PNAs bearing **F** building block (**PNA 1-3**) and **ON 2** (upper panel) and PNAs bearing **f** building block (**PNA 4-6**) and **ON 1** (lower panel). a) Activation of the furan unit; b) reaction of the keto-enal with the exocyclic amine; c) ring closure reaction; d) re-aromatization reaction, through dehydration.