

## Supporting Information

# Synthesis, Anti-proliferative Activity, and Molecular Docking Study of New Series of 1,3-5-Triazine Schiff Base Derivatives

Hessa H. Al Rasheed,\* Azizah M. Malebari, Kholood A. Dahlous, Darren Fayne, Ayman El-Faham\*

### Table of Content

Method S1	<i>In vitro</i> anti-proliferative assay
Method S2	Annexin V/PI Apoptotic Assay
Figure S1	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4a
Figure S2	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4b
Figure S3	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4c
Figure S4	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4d
Figure S5	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4e
Figure S6	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4f
Figure S7	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4g
Figure S8	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4h
Figure S9	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4i
Figure S10	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4j
Figure S11	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4k
Figure S12	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4l
Figure S13	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4m
Figure S14	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4n
Figure S15	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4o
Figure S16	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4p
Figure S17	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4q
Figure S18	<sup>1</sup> H-NMR and <sup>13</sup> C-NMR for compound 4r

### **Method S1: In vitro anti-proliferative assay**

The synthesized compounds **4a-r** were evaluated for anti-proliferative using the MTT viability assay of MCF-7 and HCT-116 cell lines and to calculate the relative IC<sub>50</sub> values for each compound). Cells were seeded in triplicate in 96-well plates at a density of  $10 \times 10^3$  cells/mL in a total volume of 200  $\mu$ L per well. 1% ethanol (v/v) was used as a vehicle control. After 24 h, they were treated with 2  $\mu$ L test compound which had been pre-prepared as stock solutions to furnish the concentration range of study, 0.1  $\mu$ M to 100  $\mu$ M, and re-incubated for a further 72 h. The culture medium was then removed and the cells washed with 100  $\mu$ L phosphate buffered saline (PBS) and 100  $\mu$ L MTT added, to reach a final concentration of 5 mg/mL MTT added. Cells were incubated for 3 h in darkness at 37 °C. At this point solubilisation was begun through the addition of 200  $\mu$ L DMSO and the cells maintained at room temperature in darkness for 20 min to ensure thorough color diffusion before reading the absorbance the optical density was detected with a microplate reader at 570 nm. Results were expressed as percentage viability relative to vehicle control (100%). Dose response curves were plotted and IC<sub>50</sub> values (concentration of drug resulting in 50% reduction in cell survival) were obtained using the commercial software package Prism (GraphPad Software, Inc., La Jolla, CA, USA). All the experiments were repeated in at least three independent experiments.

### **Method S2: Annexin V/PI Apoptotic Assay**

Apoptotic cell death was detected by flow cytometry using Annexin V and propidium iodide (PI). MCF-7 Cells were seeded in 6 well plated at density of  $1 \times 10^5$  cells/mL and treated with vehicle (0.1 % (v/v) DMSO) and compound **4b** (6  $\mu$ M) for 24 hr. Cells were then harvested and prepared for flow cytometric analysis. Cells were washed in 1X binding buffer (20X binding buffer: 0.1M HEPES, pH 7.4; 1.4 M NaCl; 25 mM CaCl<sub>2</sub> diluted in dH<sub>2</sub>O) and incubated in the dark for 30 minutes on ice in Annexin V-containing binding buffer [1:100]. Cells were then washed once in binding buffer and then re-suspended in PI-containing binding buffer [1:1000]. Samples were analysed immediately using the BD accuri flow cytometer and prism software for analysis the data. Four populations are produced during the assay Annexin V and PI negative (Q4, healthy cells), Annexin V positive and PI negative (Q3, early apoptosis), Annexin V and PI positive (Q2, late apoptosis) and Annexin V negative and PI positive (Q1, necrosis).

Figure S1:  $^1\text{H-NMR}$  and  $^{13}\text{C-NMR}$  for compound 4a





Figure S2:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound 4b





Figure S3:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound 4c





Figure S4:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound 4d





Figure S5:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound 4e





Figure S6:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound 4f





Figure S7:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound 4g




Figure S8:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound 4h





Figure S9:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound **4i**





Figure S10:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound 4j




Figure S11:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound **4k**





Figure S12:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound 4l





Figure S13:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound 4m





Figure S14:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound 4n





Figure S15:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound 4o





Figure S16:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound 4p




Figure S17:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound 4q




Figure S18:  $^1\text{H}$ -NMR and  $^{13}\text{C}$ -NMR for compound 4r

