



Article

# Strategy for Accurate Detection of *Escherichia coli* O157:H7 in Ground Pork Using a Lateral Flow Immunoassay

Song Cheng  $^1$ , Ming-hui Chen  $^1$ , Gang-gang Zhang  $^1$ , Zhi-biao Yu  $^2$ , Dao-feng Liu  $^1$ , Yong-hua Xiong  $^1$ , Hua Wei  $^1$  and Wei-hua Lai  $^{1,*}$ 

- State Key Laboratory of Food Science and Technology, Nanchang University, Nanchang 330047, China; 15797769982@163.com (S.C.); minghuiai@sina.com (M.-h.C.); joezgg@sina.com (G.-g.Z.); defoelau@163.com (D.-f.L.); xiongyonghua@ncu.edu.cn (Y.-h.X.); weihua@ncu.edu.cn (H.W.)
- Department of Science and Technology, Nanchang University, Nanchang 330047, China; yzb1156@126.com
- \* Correspondence: laiweihua@ncu.edu.cn; Tel.: +86-138-7917-8802

Academic Editors: Jean-Louis Marty, Silvana Andreescu and Akhtar Hayat Received: 9 February 2017; Accepted: 28 March 2017; Published: 2 April 2017

**Abstract:** *Escherichia coli* O157:H7 is known to cause serious diseases including hemorrhagic colitis and hemolytic uremic syndrome. A gold nanoparticle lateral flow immunoassay (Au-LFIA) was used to detect *Escherichia coli* O157:H7 in ground pork samples. False-positive results were detected using Au-LFIA; a *Citrobacter freundii* strain was isolated from the ground pork samples and identified by using CHROmagar<sup>TM</sup> plates, API 20E, and 16S RNA sequencing. Since *C. freundii* showed cross-reactivity with *E. coli* O157:H7 when Au-LFIA test strips were used, a novel method combining modified enrichment with a lateral flow immunoassay for accurate and convenient detection of *E. coli* O157:H7 in ground pork was developed in this study to minimize these false positives. MacConkey broth was optimized for *E. coli* O157:H7 enrichment and *C. freundii* inhibition by the addition of 5 mg/L potassium tellurite and 0.10 mg/L cefixime. Using the proposed modified enrichment procedure, the false-positive rate of ground pork samples spiked with 100 CFU/g *C. freundii* decreased to 5%.

Keywords: enrichment; E. coli O157:H7; ground pork sample; false positive

## 1. Introduction

Escherichia coli O157:H7 is a dangerous foodborne pathogen because of its low infectious dose (minimum 10 cells) and high pathogenicity (diarrhea, hemorrhagic colitis, and hemolytic uremic syndrome) [1]. The "gold standard" for detecting *E. coli* O157:H7 in food samples is a traditional separation-identification method which is time-consuming (at least five days) and laborious (involves pre-enrichment, selective enrichment, culture isolation, and identification). Other methods, including the polymerase chain reaction (PCR) [2,3] and enzyme-linked immunosorbent assay (ELISA) [4,5], also may require laborious procedures and expensive instruments.

Lateral flow immunoassay (LFIA), which is widely used in the field of food safety, presents a number of advantages, including ease of use, rapidity, and sensitivity [6]. However, LFIA-based methods have the disadvantage that the antibody used in the immunoassay can show cross-reactivity with bacteria other than the target [7,8].

A gold nanoparticle LFIA (Au-LFIA) method was previously optimized for detecting *Escherichia coli* O157:H7 with high sensitivity [9–15]. Detection of *E. coli* O157:H7 in uninoculated ground pork samples using Au-LFIA test strips prepared in our laboratory in combination with the enrichment in modified *E.Coli* (EC) broth yielded positive results. In this study, the reasons behind the false-positive results obtained when Au-LFIA test strips are used to detect *E. coli* O157:H7 in ground pork were determined.

Sensors **2017**, 17, 753

A novel method combining modified enrichment with Au-LFIA for accurate detection of *E. coli* O157:H7 in ground pork was developed to minimize false positives (Figure 1).




Figure 1. Overall process of the immunochromatographic assay for detecting E. coli O157:H7.

# 2. Materials and Methods

## 2.1. Materials

Bacterial strains: The five *E. coli* O157:H7 strains and 40 non–*E. coli* O157:H7 strains used in this study are described in Table 1. Twenty ground pork samples were purchased from Guohong Pork Slaughterhouse (Nanchang, China), and another 20 ground pork samples were purchased from local supermarkets in Nanchang, China. A homogenizer and stomacher bags were purchased from Voshin Instruments Co., Ltd. (Wuxi, China). Nitrocellulose membrane, sample pad, conjugate release pad, and absorbent pad were obtained from Millipore (Bedford, MA, USA). Anti-*E. coli* O157:H7 monoclonal antibody (mAb, 10C5-H3-B6), anti-*E. coli* O157:H7 polyclonal antibodies (pAb), and goat anti-mouse antibody pAb were obtained from Wuxi Zodoboer Biotech. Co., Ltd. (Wuxi, China). Novobiocin, modified EC (mEC) culture medium, potassium tellurite, cefixime, and MacConkey Broth (CT-SMAC) culture medium were purchased from Beijing Land Bridge Technology Co., Ltd. (Beijing, China). Colloidal gold test strip reader was obtained from Suzhou Helmen Precision Instrument Co., Ltd. (Suzhou, China). CHROmagar<sup>TM</sup> O157 agar for *E. coli* O157 was purchased from CHROMagar (Paris, France). API 20E test strips were purchased from Biomerieux (Lyon, France). Then 16S RNA sequencing was performed by Genscipt Company (Nanjing, China).

## 2.2. Preparation of Au-LFIA Test Strips

Au-LFIA test strips to detect *E. coli* O157:H7 were prepared in laboratory as described previously [14]. The test strip consisted of sample pad, conjugate pad, nitrocellulose membrane (NC), and absorbent pad. Anti-*E. coli* O157:H7 polyclonal antibody was applied to the nitrocellulose membrane as the test line. Goat anti-mouse antibody was applied to the nitrocellulose membrane as the control lines. Anti-*E. coli* O157:H7 monoclonal antibody-AuNPs complex was applied to a conjugate pad.

# 2.3. Pretreatment of Ground Pork Samples

Twenty ground pork samples were tested negative for  $E.\ coli\ O157:H7$ . Then 25 g of these ground pork samples was transferred into a stomacher bag as negative sample. Another 25 g of these ground pork samples was transferred into a stomacher bag and inoculated 2 CFU/g  $E.\ coli\ O157:H7$  as positive samples.

Sensors **2017**, 17, 753

# 2.4. Enrichment with Modified EC Broth and Evaluation by Using Au-LFIA Test Strip

Twenty negative samples and 20 positive samples (25 g) were mixed with 225 mL of modified EC broth containing 20 mg/L of novobiocin respectively, and stomached (Seward 400 Stomacher, Norfolk, UK) for 2 min. All of the samples were incubated at 37 °C with shaking at 160 r/min for 12 h. One hundred microliters of the 40 enriched broths of the ground pork samples was respectively added to the sample pad of Au-LFIA test strips for detection. When the target analyte (*E. coli* O157:H7) was added onto the sample pad of the test strip, it flowed laterally through the test strip. The *E. coli* O157:H7 interacted with the antibody-AuNPs complex in the conjugate pad, and aggregated subsequently at the T line because of the specific interaction between the *E. coli* O157:H7 and the antibody, which leads to a red color of the T line. Two visual bands (control line and test line) indicated a positive result. One visual band (control line) indicated a negative test result. The colored gold-antibody conjugate should bind to the control line and form a red-colored band regardless of the presence of *E. coli* O157:H7. The signal intensity of the test line can be detected by colloidal gold test strip reader. When the sample is positive, the signal intensity of the test line is equal or greater than 30.

# 2.5. Specificity of the Au-LFIA Test Strip

Five *E. coli* O157:H7 strains and 40 non- *E. coli* O157:H7 strains (Table 1) were cultured in Luria-Bertani medium (LB, Oxoid, Basingstoke, UK) at 37 °C for 20 h before use. Then 10<sup>5</sup> CFU/mL of strains were prepared by serial dilutions of cultures in phosphate buffered saline (PBS, Sigma Chemical Company, St. Louis, MO, USA, 0.01 M, pH 7.4). One hundred microliters of the 45 strains were tested by the Au-LFIA test strip for evaluating the specificity of the method. The non-*E. coli* O157:H7 strains cannot interact with the antibody-labeled AuNPs and no red line develops at the test line.

Species	Source	Result (Approx. 10 <sup>5</sup> CFU/mL)
E. coli O157:H7	ATCC 43888	+
E. coli O157:H7	CMCC 44828	+
E. coli O157:H7	NCTC 12900	+
E. coli O157:H7	XY0540 a	+
E. coli O157:H7	XY0480 a	+
Bacillus cereus	CMCC 63303	_
Bacillus cereus	CMCC 63305	_
Bacillus cereus	SLK <sup>a</sup>	_
Bacillus licheniformis	CMCC 63519	_
Bacillus subtilis	BD366 <sup>a</sup>	_
Candida albicans	ATCC 10231	_
Candida albicans	Z1 <sup>a</sup>	_
Cronobacter cloacae	CMCC 45301	_
Cronobacter sakazakii	CMCC 45401	_
Cronobacter sakazakii	CMCC 45402	_
enteropathogenic E. coli	CMCC 44496	_
E. coli	CMCC 44102	_
E. coli	ATCC 25922	_
Lactobacillus bulgaricus	F1 <sup>a</sup>	_
Listeria grayi	ATCC 25401	_
Listeria innocua	ATCC 33090	_
Listeria innocua	ATCC 11288	_
Listeria ivanovii	ATCC 19119	_
Listeria monocytogenes	ATCC 13932	_
Listeria monocytogenes	CMCC 54001	_
Listeria monocytogenes	CMCC 54007	_
Listeria welshimeri serovar 6b	ATCC 35897	

Sensors **2017**, 17, 753 4 of 8

TO 1 '	1.	•	$\sim$	
Tab:	Ie -		เดท	11

Species	Source	Result (Approx. 10 <sup>5</sup> CFU/mL)
Listeria seeligeri	ATCC 35967	_
Micrococcus luteus	CMCC 28001	_
Proteus vulgaris	CMCC 49027	_
Pseudomonas aeruginosa	ATCC 27853	_
Pseudomonas aeruginosa	CMCC 10104	_
Salmonella anatum	ATCC 9270	_
Salmonella choleraesuis	ATCC 13312	_
Salmonella choleraesuis	CICC 21493	_
Salmonella enteritidis	ATCC 13076	_
Salmonella paratyphi A	ATCC 9150	_
Salmonella enterica	ATCC 10708	_
Salmonella typhimurium	ATCC 13311	_
Serratia marcescens	CMCC 41002	_
Shigella flexneri	ATCC 12022	_
Shigella sonnei	CMCC 51592	_
Staphylococcus aureus	CMCC 26001	_
Staphylococcus aureus	CMCC 26003	_
Vibrio parahaemolyticus	CGMCC 1.1997	_

 $<sup>^{</sup>a}$  Strains were from Jiangxi Province Center for Disease Control and Prevention; "+": positive result, "-": negative result.

## 2.6. Identification of the False-Positive Bacterium

CHROmagar<sup>TM</sup> O157 agar without potassium tellurite and cefixime was used after ground pork sample enrichment. The enrichment broth was diluted 1:100,000 with sterile phosphate buffer, and the diluted broth (100  $\mu$ L) was streaked onto CHROmagar<sup>TM</sup> O157 agar. After incubation at 37 °C for 24 h, bacteria with different colors and colony morphology were inoculated into LB medium for enrichment and identified through API 20E and 16S RNA sequencing.

# 2.7. Cross Activity of C. freundii with Au-LFIA Strip

The isolated *C. freundii* strain was identified, inoculated into 250 mL of modified EC broth (10 CFU/mL), and incubated at 37  $^{\circ}$ C with shaking at 160 r/min for 12 h. The sample of the enrichment broth (100  $\mu$ L) was added to the sample pad of Au-LFIA test strips. Three experiments were repeated.

# 2.8. Optimization of Enrichment Conditions

CT-SMAC culture medium with potassium tellurite and cefixime was optimized by adjusting the concentrations of potassium tellurite (2.50 mg/L, 3.75 mg/L, and 5.00 mg/L) and cefixime (0.05 mg/L, 0.075 mg/L, and 0.10 mg/L). Two portions of the ground pork samples that had been tested to be negative for *E. coli* O157:H7 and *C. freundii* were spiked with *E. coli* O157:H7 (2 CFU/g) and *C. freundii* which was identified in 2.6 (100 CFU/g). Two hundred and twenty-five milliliters of CT-SMAC with potassium tellurite and cefixime was mixed with 25 g of the ground pork samples and stomached for 2 min. All of the test samples were incubated at 37  $^{\circ}$ C with shaking at 160 r/min for 12 h.

# 2.9. Evaluation of the Optimized Enrichment by Using Au-LFIA Test Strip

Twenty negative ground pork samples were spiked with 100 CFU/g *C. freundii* as negative samples, while another 20 negative ground pork samples were spiked with 2 CFU/g *E. coli* O157:H7 as positive samples. Twenty-five grams of these 40 samples were respectively mixed with 225 mL of the optimized CT-SMAC culture medium and stomached for 2 min. All of the samples were incubated at  $37 \,^{\circ}\text{C}$  with shaking at  $160 \,\text{r/min}$  for  $12 \,\text{h}$ . One hundred microliters of the  $40 \,\text{enrichment}$  broth of the ground pork samples were then respectively added to the sample pad of Au-LFIA test strips for detection.

Sensors **2017**, *17*, 753 5 of 8

#### 3. Results and Discussion

# 3.1. Specificity of the Au-LFIA

The results showed that the five strains of *E. coli* O157:H7 were detected successfully, and no cross-reaction with the other 40 strains was observed (Table 1).

## 3.2. False-Positive Result of Au-LFIA Test Strip

Forty inoculated positive ground pork samples (20 from the slaughterhouse and 20 from the supermarket) and 40 negative control ground pork samples (20 from the slaughterhouse and 20 from the supermarket) were enriched by using modified EC broth. The samples were then tested for the presence of *E. coli* O157:H7 by using Au-LFIA test strips. Forty samples from the slaughterhouse, including 20 positive samples and 20 negative samples, yielded positive results. Twenty positive samples from supermarket showed positive test results. As well, eight out of the 20 negative samples from the supermarket showed positive test results, whereas 12 out of 20 uninoculated samples showed negative test results (Table 2).

The results obtained from Au-LFIA test strips indicated that some bacterial species in the ground pork samples exhibit cross-reactivity with *E. coli* O157:H7 (Figure S1). All 20 negative samples from the slaughterhouse showed false positives, whereas only eight of the 20 negative samples from the supermarket revealed false positives.

Samples		Negative Result of Au-LFIA	Positive Result of Au-LFIA
Campalas from the slavebook augo	10 positive samples	0	20
Samples from the slaughterhouse	10 negative samples	0	20
Company to a second sec	10 positive samples	0	20
Samples from the supermarket	10 negative samples	12	8

**Table 2.** Results of ground pork samples detected using Au-LFIA test strips.

Positive samples are inoculated with 2 CFU/g E. coli O157:H7; negative samples are without inoculated with E. coli O157:H7.

# 3.3. Isolation and Identification of C. freundii from Incubated Broth

The diluted modified EC broth (100 µL) was streaked onto CHROmagar<sup>TM</sup> O157 agar. After incubation at 37 °C for 24 h, colonies with different colors and colony morphologies were produced. *E. coli* O157:H7 appeared mauve in CHROmagar<sup>TM</sup> plates (Figure S2A), while some non-O157 bacteria showed other colors. The bacteria in 14 negative samples with false-positive Au-LFIA results were isolated using the CHROmagar<sup>TM</sup> plate. These bacteria exhibited cross-reactivity with *E. coli* O157:H7 when the Au-LFIA test strips were used for detection and appeared blue on the CHROmagar<sup>TM</sup> plate (Figure S2B). The bacteria producing false positives were identified by API 20E (Figure S3), and one bacterium, namely *C. freundii*, was identified through 16S RNA sequencing (Supplementary Material 4). The nucleotide sequence of *C. freundii* was identified with 99.9% accuracy by NCBI Blast, and the results obtained agreed with the API 20E findings.

# 3.4. Cross-Reactivity of C. freundii with Au-LFIA Strip


The positive results indicated that the Au-LFIA test strips for detecting *E. coli* O157:H7 exhibited cross-reactivity with *C. freundii*.

# 3.5. Optimization of the Modified Culture Medium for C. freundii Inhibition

A series of various potassium tellurite and cefixime concentrations in CT-SMAC were studied to demonstrate the inhibition of *C. freundii* (Figure 2). When the potassium tellurite and cefixime

Sensors **2017**, 17, 753 6 of 8

concentrations in CT-SMAC were 5 and 0.10 mg/L, respectively, *C. freundii* did not multiply and *E. coli* O157:H7 was enriched. One hundred microliters of the enrichment broth of the ground pork samples was added to the sample pad of Au-LFIA test strips prepared in our laboratory. The signal intensities of the test lines of the strips with ground pork broth samples spiked with *E. coli* O157:H7 and *C. freundii* were 175 and 0, respectively.



Concentration of Potassium tellurite/Cefixime in CT-SMS (mg/L)

**Figure 2.** Signal intensities of the Au-LFIA strip test lines with pork broth samples spiked with  $E.\ coli$  O157:H7 and  $C.\ freundii\ (n = 3)$ .

## 3.6. Evaluation of Modified CT-SMAC in Ground Pork Test

Twenty ground pork samples that were negative for *E. coli* O157:H7 and *C. freundii* were spiked with 100 CFU/g *C. freundii* as negative controls and 2 CFU/g *E. coli* O157:H7 as positive controls. After enrichment by the developed culture medium, 100  $\mu$ L of the enrichment broth of the ground pork samples was detected using Au-LFIA test strips, as well as CHROmagar<sup>TM</sup> plates. Results (Table 3) indicated that the modified CT-SMAC was suitable for *E. coli* O157:H7 enrichment in ground pork samples. Using this enrichment procedure, ground pork samples spiked with 2 CFU/g *E. coli* O157:H7 showed 100% positive results and only a 5% false-positive result from *C. freundii*.

Previous studies had shown that some strains of *C. freundii* presented cross-reactivity with anti-O157 sera [16,17]. With the developed enrichment procedure, Au-LFIA had good specificity. Bennett and Zadik have also obtained a good specificity result with potassium tellurite and cefixime [18,19]. Some scholars also acquired sensitive and specific results based on lectin recognition of *E. coli* O157:H7 [20].

Pork Samples	CHROmagar <sup>TM</sup> Plates		Au-LFIA Test Strip	
•	+	_	+	_
Positive samples ( <i>E. coli</i> O157:H7) ( $n = 20$ )	20	0	20	0
Negative samples ( <i>C. freundii</i> ) ( $n = 20$ )	0	20	1	19

**Table 3.** Evaluation of modified CT-SMAC in ground pork samples.

## 4. Conclusions

In this study, *E. coli* O157:H7 in ground pork samples was detected by using Au-LFIA test strips. A large number of false-positive results were obtained. The *C. freundii* strain was isolated and identified from the ground pork samples and determined to induce these false positives. A modified enrichment procedure by the addition of 5 mg/L potassium tellurite and 0.10 mg/L cefixime was evaluated for the enrichment of *E. coli* O157:H7 and the inhibition of *C. freundii*. Combining the modified enrichment

<sup>&</sup>quot;+": positive result, "-": negative result.

Sensors **2017**, 17, 753 7 of 8

procedure with Au-LFIA, ground pork samples spiked with 2 CFU/g *E. coli* O157:H7 showed 100% positive results and a 5% false-positive result from *C. freundii*.

**Supplementary Materials:** The following are available online at http://www.mdpi.com/1424-8220/17/4/753/s1, Figure S1: Detection results of the Au-LFIA test strips, Figure S2: CHROmagarTM plates used to isolate bacteria from the enrichment broth. (A) E. coli O157:H7. (B) Cross-reactive bacterium, Figure S3: Identification results of the cross-reactive bacterium by API 20E, Supplementary materials 4: Nucleotide sequence of 16S RNA of Citrobacter freundii.

**Acknowledgments:** We are grateful to the Open Project Program of the State Key Laboratory of Food Science and Technology of Nanchang University (SKLF-KF-201616), Free explore issue of State Key Laboratory of Food Science and Technology of Nanchang University (SKLF-ZZA-201612), and the Jiangxi Agriculture Research System Project (JXARS-03).

**Author Contributions:** Wei-hua Lai conceived and designed the experiments; Song Cheng and Ming-hui Chen performed the experiments; Hua Wei and Yong-Hua Xiong analyzed the data; Zhi-biao Yu, Dao-feng Liu, and Gang-gang Zhang contributed reagents/materials/analysis tools; Wei-hua Lai and Song Cheng wrote the paper.

Conflicts of Interest: The authors declare no conflict of interest.

#### References

- 1. Pennington, H. Escherichia coli O157. Lancet 2010, 376, 1428–1435. [CrossRef]
- 2. Liu, Y.; Mustapha, A. Detection of viable *Escherichia coli* O157:H7 in ground beef by propidium monoazide real-time PCR. *Int. J. Food Microbial.* **2014**, 170, 48–54. [CrossRef] [PubMed]
- 3. Peng, T.; Zhang, F.S.; Yang, W.C.; Li, D.X.; Chen, Y.; Xiong, Y.H.; Wei, H.; Lai, W.H. Lateral-Flow Assay for Rapid Quantitative Detection of Clorprenaline Residue in Swine Urine. *J. Food Prot.* **2014**, 77, 1824–1829. [CrossRef] [PubMed]
- 4. Amani, J.; Ahmadpour, A.; Fooladi, A.A.I.; Nazarian, S. Detection of *E. coli* O157:H7 and *Shigella dysenteriae* toxins in clinical samples by PCR-ELISA. *Braz. J. Infect. Dis.* **2015**, *19*, 278–284. [CrossRef] [PubMed]
- 5. Deng, S.L.; Shan, S.; Xu, C.L.; Liu, D.F.; Xiong, Y.H.; Wei, H.; Lai, W.H. Sample Preincubation Strategy for Sensitive and Quantitative Detection of Clenbuterol in Swine Urine Using a Fluorescent Microsphere–Based Immunochromatographic Assay. *J. Food Prot.* **2014**, 77, 1998–2003. [CrossRef] [PubMed]
- 6. Shan, S.; Lai, W.; Xiong, Y.; Wei, H.; Xu, H. Novel strategies to enhance lateral flow immunoassay sensitivity for detecting foodborne pathogens. *J. Agric. Food Chem.* **2015**, *63*, 745–753. [CrossRef] [PubMed]
- 7. Borrebaeck, C.A. Antibodies in diagnostics–from immunoassays to protein chips. *Immunol. Today* **2000**, *21*, 379–382. [CrossRef]
- 8. Van Der Veer, J.; Lewis, R.J.; Emtiazjoo, A.M.; Allen, S.D.; Wheat, L.J.; Hage, C.A. Cross-reactivity in the Platelia™ Aspergillus enzyme immunoassay caused by blastomycosis. *Med. Mycol.* **2012**, *50*, 396–398. [CrossRef] [PubMed]
- 9. Song, C.; Liu, C.; Wu, S.; Li, H.; Guo, H.; Yang, B.; Qiu, S.; Li, J.; Liu, L.; Zeng, H. Development of a lateral flow colloidal gold immunoassay strip for the simultaneous detection of *Shigella boydii* and *Escherichia coli* O157: H7 in bread, milk and jelly samples. *Food Control* **2016**, *59*, 345–351. [CrossRef]
- Cho, I.-H.; Bhunia, A.; Irudayaraj, J. Rapid pathogen detection by lateral-flow immunochromatographic assay with gold nanoparticle-assisted enzyme signal amplification. *Int. J. Food Microbiol.* 2015, 206, 60–66. [CrossRef] [PubMed]
- 11. Xie, Q.-Y.; Wu, Y.-H.; Xiong, Q.-R.; Xu, H.-Y.; Xiong, Y.-H.; Liu, K.; Jin, Y.; Lai, W.-H. Advantages of fluorescent microspheres compared with colloidal gold as a label in immunochromatographic lateral flow assays. *Biosens. Bioelectron.* **2014**, *54*, 262–265. [CrossRef] [PubMed]
- 12. Cui, X.; Xiong, Q.-R.; Xiong, Y.-H.; Shan, S.; Lai, W.H. Establishing of a method combined immunomagnetic separation with colloidal gold lateral flow assay and its application in rapid detection of *Escherichia coli* O157:H7. *Chin. J. Anal. Chem.* **2013**, *41*, 1812–1816. [CrossRef]
- 13. Wang, J.-Y.; Chen, M.-H.; Sheng, Z.-C.; Liu, D.-F.; Wu, S.-S.; Lai, W.-H. Development of colloidal gold immunochromatographic signal-amplifying system for ultrasensitive detection of *Escherichia coli* O157:H7 in milk. *RSC Adv.* **2015**, *5*, 62300–62305. [CrossRef]
- 14. Cui, X.; Huang, Y.; Wang, J.; Zhang, L.; Rong, Y.; Lai, W.; Chen, T. A remarkable sensitivity enhancement in a gold nanoparticle-based lateral flow immunoassay for the detection of *Escherichia coli* O157:H7. *RSC Adv.* **2015**, *5*, 45092–45097. [CrossRef]

Sensors **2017**, 17, 753

15. Chen, M.; Yu, Z.; Liu, D.; Peng, T.; Liu, K.; Wang, S.; Xiong, Y.; Wei, H.; Xu, H.; Lai, W. Dual gold nanoparticle lateflow immunoassay for sensitive detection of *Escherichia coli* O157:H7. *Anal. Chim. Acta* **2015**, *876*, 71–76. [CrossRef] [PubMed]

- 16. Sowers, E.G.; Wells, J.G.; Strockbine, N.A. Evaluation of commercial latex reagents for identification of O157 and H7 antigens of *Escherichia coli. J. Clin. Microbiol.* **1996**, *34*, 1286–1289. [PubMed]
- 17. Navarro, A.; Eslava, C.; de la Torre, G.G.; Leon, L.A.; Licona, D.; Leon, L.; Zarco, L.A.; Cravioto, A. Common epitopes in LPS of different Enterobacteriaceae are associated with an immune response against *Escherichia coli* O157 in bovine serum samples. *J. Med. Microbial.* **2007**, *56*, 1447–1454. [CrossRef] [PubMed]
- 18. Bennett, A.R.; Macphee, S.; Betts, R.P. Evaluation of methods for the isolation and detection of *Escherichia coli* O157 in minced beef. *Lett. Appl. Microbiol.* **1995**, 20, 375. [CrossRef] [PubMed]
- 19. Zadik, P.M.; Chapman, P.A.; Siddons, C.A. Use of tellurite for the selection of verocytotoxigenic *Escherichia coli* O157. *J. Med. Microbiol.* **1993**, 39, 155–158. [CrossRef] [PubMed]
- 20. Wang, Y.; Ye, Z.; Si, C.; Ying, Y. Monitoring of *Escherichia coli* O157:H7 in food samples using lectin based surface plasmon resonance biosensor. *Food Chem.* **2013**, *136*, 1303–1308. [CrossRef] [PubMed]



© 2017 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (http://creativecommons.org/licenses/by/4.0/).