

Supplementary Materials: Matrix-Matched Iron-Oxide Laser Ablation ICP-MS U-Pb Geochronology Using Mixed Solution Standards

Liam Courtney-Davies, Zhiyong Zhu, Cristiana L. Ciobanu, Benjamin P. Wade,
Nigel J. Cook, Kathy Ehrig, Alexandre R. Cabral and Allen Kennedy

Supplementary Text File S1. Additional analytical procedures for calculating and certifying the Pb/U ratios.

The mixed solution standard was made from single natural U and Pb solutions. The concentrations of these two solutions were calibrated using Q-ICP-MS at Adelaide Microscopy.

- Q-ICP-MS results: Pb concentration of the mother bottle = 9985 ppb; U concentration of the mother bottle = 9994 ppb.

The Pb isotope ratios were analyzed by MC-ICP-MS at Nanjing University.

- MC-ICP-MS results: $^{206}Pb/^{204}Pb = 17.9895 \pm 0.0005$; $^{207}Pb/^{204}Pb = 15.5847 \pm 0.0005$; $^{208}Pb/^{204}Pb = 38.0495 \pm 0.0013$; $^{207}Pb/^{206}Pb = 0.86632 \pm 0.00004$.

The constant $^{238}U/^{235}U$ (137.88), in combination with the Pb isotope ratio results, allows us to calculate $^{206}Pb/^{238}U$.

To achieve this, we must calculate the atomic weight of the Pb and U solution standards:

X = Proportion of isotopes; M = Mass of isotope

- Pb atomic weight = $X^{204} \times M^{204} + X^{206} \times M^{206} + X^{207} \times M^{207} + X^{208} \times M^{208}$
- U atomic weight = $X^{235} \times M^{235} + X^{238} \times M^{238}$

To then calculate the ratio of ^{206}Pb and ^{238}U , the correct amount of Pb and U standard solution to create a concentration compatible with the sample to be analyzed must be weighed.

Con = Concentration of element; W = Weight of solution

$$\frac{Pb_{206}}{U_{238}} = \frac{\frac{Con_{pb} \times W_{Pb-sol}}{Pb_{atomic\ weight}} \times X_{206}}{\frac{Con_U \times W_U}{U_{atomic\ weight}} \times X_{238}} \quad (1)$$

$$\frac{Pb_{207}}{U_{235}} = \frac{\frac{Con_{pb} \times W_{Pb-sol}}{Pb_{atomic\ weight}} \times X_{207}}{\frac{Con_U \times W_U}{U_{atomic\ weight}} \times X_{235}} \quad (2)$$

$^{207}Pb/^{206}Pb = 0.86632 \pm 0.00004$ (measured by MC-ICP-MS)

$^{206}Pb/^{238}U = 0.286356$

The following tables contain the Pb isotope ratios and U-Pb concentrations.

Table S1. Pb isotope ratios of Pb standard solution.

Name	Pb ^{206/204}	Pb ^{206/204} SE (abs)	Pb ^{207/204}	Pb ^{207/204} SE (abs)	Pb ^{208/204}	Pb ^{208/204} SE (abs)
NIST987-01	16.939	0.001	15.488	0.001	36.693	0.002
NIST987-02	16.937	0.001	15.487	0.001	36.690	0.002
Pb-1	17.990	0.001	15.585	0.001	38.049	0.001
Pb-2	17.991	0.001	15.586	0.001	38.052	0.001
Pb-3	17.988	0.000	15.583	0.000	38.047	0.001
NIST987-03	16.938	0.001	15.488	0.001	36.692	0.002
Average of Samples	17.990	0.001	15.585	0.001	38.050	0.001

Table S2. U and Pb concentration of standard solutions. Conc. = Concentration.

Type	Sample Name	204Pb		206Pb		207Pb		208Pb		235U		238U		115 In (ISTD)			Average Pb		Average U		Pb/U
		Conc.	RSD	Conc.	RSD	Conc.	RSD	Conc.	RSD	Conc.	RSD	Conc.	RSD	Conc.	Recovery %	Conc.	RSD	Conc.	RSD	Conc.	RSD
CalBlk	Blank	0.00	N/A	0.00	N/A	0.00	N/A	0.00	N/A	0.00	N/A	-	-	100	-	-	-	-	-	-	
CalStd	Standard-1 ppb	0.98	8.07	0.95	2.34	0.96	4.70	0.95	2.01	0.98	17.39	1.00	1.95	-	-	99	0.95	0.19	0.99	1.84	0.96
CalStd	Standard-2 ppb	1.97	6.05	1.97	1.74	1.96	3.03	2.00	1.88	2.04	6.45	2.02	0.84	-	-	97	1.98	1.14	2.03	0.40	0.97
CalStd	Standard-5 ppb	5.06	3.76	5.01	1.92	4.98	2.46	5.07	2.24	4.95	5.22	5.11	2.66	-	-	100	5.02	0.86	5.03	2.31	1.00
CalStd	Standard-10 ppb	9.89	4.49	9.91	3.38	9.90	3.82	10.02	3.97	9.85	4.14	9.99	3.18	-	-	103	9.94	0.65	9.92	1.02	1.00
CalStd	Standard-20 ppb	19.42	2.66	19.63	2.47	19.58	2.62	19.98	2.17	19.34	3.52	19.97	1.80	-	-	106	19.73	1.09	19.65	2.29	1.00
CalStd	Standard-50 ppb	50.25	1.82	50.17	2.09	50.19	2.63	51.55	2.91	50.30	2.23	51.82	2.48	-	-	104	50.64	1.56	51.06	2.11	0.99
Sample	Blank	<0.000	N/A	<0.000	N/A	<0.000	N/A	<0.000	N/A	<0.000	N/A	<0.000	N/A	-	-	103	-	-	-	-	-
Sample	1	1.17	12.14	1.04	5.38	1.03	4.95	1.04	5.54	1.01	14.29	0.99	5.38	-	-	105	1.03	0.78	1.00	1.98	1.04
Sample	2	4.50	4.97	4.07	2.91	3.95	2.33	4.11	2.89	3.64	4.42	3.77	1.29	-	-	115	4.04	5.95	3.71	2.50	1.09
Sample	3	11.80	2.90	11.88	2.78	11.79	1.75	12.05	2.11	11.37	4.19	11.76	2.39	-	-	111	11.91	1.09	11.56	2.41	1.03
Sample	4	43.61	1.65	44.42	2.10	44.09	2.21	44.97	2.38	38.70	3.25	40.39	2.95	-	-	108	44.50	1.00	39.54	3.01	1.13

Table S3. Relation coefficient of standards.

Tune Step	Mass	Name	Relation Coefficient	Detection Limit
2	204	Pb	0.99990	0.10246
2	206	Pb	0.99996	0.06304
2	207	Pb	0.99995	0.07154
2	208	Pb	0.99999	0.06215
2	235	U	0.99987	0.09124
2	238	U	0.99998	0.03269



© 2016 by the authors; licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons by Attribution (CC-BY) license (<http://creativecommons.org/licenses/by/4.0/>).