

Supplementary Data

Table S1. Targets of the Novodiag® Stool Parasites (NSP) assay.

Protozoa targeted by the NSP
<i>Balantidium coli</i>
<i>Blastocystis</i> spp.
<i>Cryptosporidium</i> spp. ^a
<i>Cyclospora cayetanensis</i>
<i>Cystoisospora belli</i>
<i>Dientamoeba fragilis</i>
<i>Entamoeba histolytica</i>
<i>Giardia intestinalis</i>
Microsporidia targeted by the NSP
<i>Encephalitozoon</i> spp. ^b
<i>Enterocytozoon bieneusi</i>
Helminths targeted by the NSP
<i>Ancylostoma duodenale</i>
<i>Ascaris lumbricoides</i> / suum
<i>Clonorchis sinensis</i> / <i>Opisthorchis</i> spp. / <i>Metorchis</i> spp.
<i>Diphyllobothrium latum</i> / <i>nihonkaiense</i>
<i>Enterobius vermicularis</i>
<i>Fasciola</i> spp. ^c
<i>Fasciolopsis buski</i>
<i>Hymenolepis diminuta</i>
<i>Hymenolepis nana</i>
<i>Necator americanus</i>
<i>Schistosoma mansoni</i>
<i>Schistosoma</i> spp. ^d
<i>Strongyloides stercoralis</i>
<i>Taenia saginata</i> / <i>asiatica</i>
<i>Taenia solium</i>
<i>Trichuris</i> spp.

^a — Detection of at least *C. hominis*, *C. meleagridis*, *C. parvum*, *C. ubiquitum* and *Cryptosporidium* spp. Chipmunk genotype I. ^b — Detection of at least *E. cuniculi*, *E. hellem* and *E. intestinalis*. ^c — Detection of *F. hepatica* and *F. gigantica*. ^d — Detection of at least *S. bovis*, *S. curassoni*, *S. edwardiense*, *S. haematobium*, *S. hippopotami*, *S. incognitum*, *S. intercalatum*, *S. japonicum*, *S. leiperi*, *S. malayensis*, *S. mansoni*, *S. margrebowiei*, *S. mattheei*, *S. mekongi* and *S. sinensium*.