



Article

Chiral recognition by DNA-immobilized TLC plate

Masanori Yamada* and Mami Inoue

Department of Chemistry, Faculty of Science, Okayama University of Science, Ridaicho 1-1, Kita-ku, Okayama 700-0005, Japan; inoue_ous@yahoo.co.jp

* Correspondence: myamada@chem.ous.ac.jp; Tel.: +81-86-286-9550




Figure S1.Photograph of L- and D-form amino acids (Asp, Val, and Arg) on DNA-TLC plate. The mobile phase is entry 1 in Table 1. The amino acid was detected by the ninhydrin reaction.




Figure S2. *R*_f values of L- and D-form amino acids on DNA-TLC plate. The mobile phase is entry 1 in Table 1. The amino acid was detected by the ninhydrin reaction. Each of the values represents the mean of three separate determinations.