

**Table S1.** Supplements used to prepare DMEM and Green's cell culture media. Reagents were purchased from Sigma (UK) unless specified.

Supplement	DMEM (6546)	Green's media
Foetal calf serum (Biowest Biosera)		10% v/v
Penicillin- streptomycin	100 IU/mL - 100 mg/mL	
Amphotericin B	0.625 µg/mL	
Ham's F12 (N4888)		3:1 v/v
DMEM AQ media (D0819)	N/A	DMEM AQ : Ham's F12
L-glutamine	2 mM	N/A
Bovine insulin		5 µg/mL
Adenine		6.25 µg/mL
T/T		10 ng/mg
Hydrocortisone	N/A	0.4 µg/mL
EGF		10 ng/mL
Cholera toxin		8.5 ng/mL