

Article

# A High-Throughput UHPLC-QqQ-MS Method for Polyphenol Profiling in Rosé Wines

Marine Lambert <sup>1</sup>, Emmanuelle Meudec <sup>1</sup>, Arnaud Verbaere <sup>1</sup>, Gérard Mazerolles <sup>1</sup>, Jérémie Wirth <sup>1</sup>, Gilles Masson <sup>2</sup>, Véronique Cheynier <sup>1</sup> and Nicolas Sommerer <sup>1,\*</sup>

- <sup>1</sup> INRA, UMR 1083 Sciences Pour l'Oenologie, Polyphenol Analysis Facility, 2 place Viala, F-34060 Montpellier, France; E-Mails: eniram.lambert@laposte.net (M.L.); meudec@supagro.inra.fr (E.M.); verbaere@supagro.inra.fr (A.V.); pfp@supagro.inra.fr (G.M.); wirth.clap@wanadoo.fr (J.W.); cheynier@supagro.inra.fr (V.C.)
- <sup>2</sup> Centre de Recherche et d'Expérimentation sur le Vin Rosé, 70 avenue Wilson, F-83550 Vidauban, France; E-Mail: gilles.masson@vignevin.com
- \* Author to whom correspondence should be addressed; E-Mail: sommerer@supagro.inra.fr; Tel.: +33-499-612-584; Fax: +33-499-612-857.

Academic Editors: Susan E. Ebeler and Helene Hopfer

Received: 13 March 2015 / Accepted: 21 April 2015 / Published: 30 April 2015

**Abstract:** A rapid, sensitive and selective analysis method using Ultra High Performance Liquid Chromatography coupled to triple-quadrupole Mass Spectrometry (UHPLC-QqQ-MS) has been developed for the quantification of polyphenols in rosé wines. The compound detection being based on specific MS transitions in Multiple Reaction Monitoring (MRM) mode, the present method allows the selective quantification of up to 152 phenolic and two additional non-phenolic wine compounds in 30 min without sample purification or pre-concentration, even at low concentration levels. This method was repeatably applied to a set of 12 rosé wines and thus proved to be suitable for high-throughput and large-scale metabolomics studies.

Keywords: rosé wine; polyphenols; UHPLC; Multiple Reaction Monitoring—MRM

### 1. Introduction

Rosé is a highly diversified type of wine widely produced and consumed worldwide. Polyphenols are key components of wines, as they are responsible for colour, taste and ultimately the quality of the wines. They constitute a widespread class of plant metabolites which can be found in grapes and can be divided into several families, for instance benzoic acids, hydroxycinnamic acids, stilbenes, flavonols, flavan-3-ols, dihydroflavonols and anthocyanins. Their occurrence in wine depends on several factors such as grape variety, wine-making process, geographic origin and age of the wine. Contrary to red wine, rosé wine polyphenolic composition has not been widely investigated and only a few studies have been published so far [1–3]. (Ultra) High Performance Liquid Chromatography ((U)HPLC) coupled to Diode Array Detection (DAD) and/or Mass Spectrometry (MS) detection has been widely used for polyphenol determination in all kinds of matrices (wine [1–24], beer [25], fruits and vegetables [22,26,27], tea [14,21,28,29], cocoa [30,31], olive oil [32], plasma [33], urine [34]). Analytical techniques used for food and/or wine polyphenolic analysis have been reviewed [7,35]. The authors detailed several separation techniques used for the fractionation and identification of each class of food polyphenols such as HPLC with UV-Visible or fluorescence detection, often hyphenated to MS detection, (High-Speed) Counter-Current Chromatography ((HS)CCC), Nuclear Magnetic Resonance (NMR) and capillary electrophoresis, sometimes coupled to MS detectors. In both reviews, the authors concluded that LC-MS methods are best suited to identification and routine quantification analyses, tandem MS being highly relevant for structural elucidation and quantification of low-level constituents.

The first polyphenol composition studies resorted to HPLC-DAD methods, for which compound quantification was based on peak areas at selected wavelengths. For instance, Dugo *et al.* [8] reported the quantification of 13 phenolic compounds by two-dimensional HPLC-DAD in 60 min. Wirth *et al.* [3] used a 55-min gradient elution program for the quantitative analysis of 28 polyphenols in rosé wine. Up to 48 phenolic compounds were quantified in red wine by Gomez-Alonso *et al.* [13], however they had to resort to both UV and fluorescence detections and to a ternary solvent chromatographic gradient of 86 min.

Analysis times could be shortened, at the expense of the diversity of the quantified polyphenols. Puértolas *et al.* [2] focused on the quantitative analysis of 14 anthocyanins in rosé wines by HPLC-DAD in a 100-min run time. *trans*-Resveratrol concentrations were determined by Bravo *et al.* [4] by HPLC-DAD (also coupled to fluorescence and electrochemical detections) in a 30-min chromatographic separation following prior sample extraction and concentration.

The use of mass spectrometry coupled to (U)HPLC-DAD techniques has brought significant improvements in polyphenol quantification, mainly by reducing analysis times and increasing sensitivity. However, most methods applied to polyphenol sub-classes analysed separately. The use of triple quadrupole-mass spectrometry (QqQ-MS) detection was relevantly described for the rapid analysis of procyanidins and anthocyanins [22,33]. Ortega *et al.* [30] showed that the use of UHPLC rather than HPLC allowed them to reduce the analysis time from 50 min to 12.5 min for cocoa samples, a less complex matrix mainly composed of procyanidins and alkaloids. Stilbene quantification has been improved regarding the analysis time and the number of quantified analytes [5,25]. Hydroxycinnamic acids were also quantified in red wines by Buiarelli *et al.* [6] by HPLC-QqQ-MS with still a rather long chromatographic separation of 70 min. Gruz *et al.* [14] reported the analysis of 17 phenolic acids by

UHPLC-QqQ-MS in 11 min in beverages *i.e.*, white wines, grapefruit juices and green tea infusions. No sample pre-concentration was required for this method, but different purification steps were necessary. A comparative study performed by Spáčil *et al.* [21] on the quantitative analysis in grape, wines and tea of 34 phenolic substances divided into four categories highlights the efficiency of the use of QqQ-MS and UHPLC. Nevertheless, these four classes of compounds were analysed separately. HPLC-DAD coupled to an Ion Trap-Mass Spectrometer (IT-MS) was also used [36] for flavonoid identification, but analysis time remained long considering the number of compounds analysed (11 flavonoid aglycones detected in 70 min).

Still, some more extensive and useful polyphenolic analysis studies have been carried out: Jaitz *et al.* [16] quantified 11 polyphenolic compounds from red wines in less than 13 min by HPLC-QqQ-MS. Twenty eight polyphenols (four procyanidins, 10 phenolic acids and 14 flavones) were quantified in cocoa extracts by HPLC-QqQ-MS in 20 min (Ortega *et al.* [31]). Vaclavik *et al.* [23] resorted to untargeted UHPLC-Q-TOF-MS for establishing a classification of red wines. Although 25 markers were identified in a 24-min gradient time run twice (in positive and negative modes), these compounds were not specifically detected nor quantified. The Mattivi group developed a method focusing on 92 anthocyanins and pigments in red wine after two preliminary semi-preparative chromatography steps (Arapitsas *et al.* [37]) that was recently applied to grape (C. Ehrhardt *et al.* [38]). This group also developed an untargeted method that allowed the detection of 57 polyphenols (among 132 searched polyphenols) in filtered red wine with a 15 min gradient UHPLC-QqQ-MS analysis (U. Vrhovsek *et al.* [39]).

Aware of all these advances in polyphenol composition analyses, a review about analytical techniques applied to wine analysis [7] points out the wine polyphenols analysis requires several specific methods due to the complexity of wine phenolics. The present work challenges this conclusion by describing a high-throughput UHPLC-QqQ-MS quantification method for rosé wine polyphenol analysis. This method was based on Multiple Reaction Monitoring (MRM) detection mode in order to increase both selectivity and sensitivity while reducing analysis times. The suitability of this method for large-scale wine analyses and metabolomics studies was demonstrated on a set of 12 rosé wines.

## 2. Results and Discussion

#### 2.1. Choice of the MRM Transitions

Twelve rosé wines were selected according to their geographic origin, their colour and their grape varieties so as to obtain a set as diverse as possible. These wines were first analysed by untargeted UHPLC-Ion Trap (IT)-MS for compound identification.

Low-concentrated compounds which were not detected by UHPLC-IT-MS were specifically targeted by UHPLC-QqQ-MRM based on transitions reported in literature. This was the case for anthocyanin ethyl-bridged proanthocyanidins [27], flavanol-anthocyanin adducts [3,20], caftaric-anthocyanins adducts [3] and anthocyanin derived pigments *i.e.*, pyranoanthocyanins [3,9,10,15,40–42] and their flavanol adducts [43]. Once the compounds identified, MRM transitions parameters for each compound were optimised on the UHPLC-QqQ-MS system and compared to data reported in the literature [5,6,14,16,17,22,25,26,28–35,44–46].

The MRM transitions were chosen either for a neutral fragment loss or an ion fragment generic to a class of polyphenols, the compound specificity being given by the parent ion mass. Each compound was detected with one quantification transition (quantifier ion) and, when available, one or two confirmation transitions (qualifier ions, see Table 1) to assess detection and quantification specificities. Moreover, the quantification transition has to produce an intense quantifier ion, so as to obtain high signal-to-noise ratios and thus reach low LOQ. When reference standards were available, the MRM transition parameters were optimised by using the Intellistant tool of the MassLynx software which consisted in automatically detecting the major fragments and optimizing cone voltages and collision energies. For other compounds, these parameters were optimised directly on wine, grape juice and Oligomeric ProCyanidins (OPC) powder samples or on enriched extracts (for anthocyanin 3-O-caffeoyl-glycosides). The optimised MRM transition parameters (ion mode, precursor and product ions, retention times) are summarized in Table 1.

Acids were characterized by the loss of the carboxylic acid group (-COO; -44 Th). For acids with a methoxy moiety, such as syringic acid and ferulic acid, the COO loss was also accompanied (or replaced) by the loss of the methyl group (-CH<sub>3</sub>; -15 Th; -COO and -CH<sub>3</sub>; -59 Th). Vanillic acid displayed a -76 Th loss corresponding to the carboxylic acid group and the -OCH<sub>3</sub> moiety. The loss of the carboxylic acid group together with a 28 a.m.u. group (probably a -CO or an ethyl moiety formed after a rearrangement) was additionally observed for protocatechuic acid and p-coumaric acid. Gallic acid also lost H<sub>2</sub>O in addition to the loss of the carboxylic acid group and of the 28 a.m.u. moiety. Grape Reaction Product (GRP), caftaric, coutaric and fertaric acids were distinguished by their tartrate fragment (149 Th), the latter three ones being characterized by their fragmentation into caffeic, p-coumaric and ferulic acids respectively (which could further undergo -COO and -CH<sub>3</sub> losses). Ethylic esters of benzoic and hydroxycinnamic acids displayed a -28 a.m.u. loss corresponding to an ethylene group, and a -73 a.m.u. loss representing the -CO-OCH<sub>2</sub>CH<sub>3</sub> group.

Cis- and trans-resveratrol were characterized by the loss of a phenol group (-94 a.m.u.), and by the formation of a hydroxytropylium ion (m/z 107 Th), and so was viniferin. *Cis-* and *trans-*piceid lost their glucose entity to form *cis-* and *trans-*resveratrol, which underwent further fragmentation as described.

Non-aglycone flavonols detection was based on the loss of their glucoside (-162 Th) or glucuronide (-176 Th) group. The main fragment (m/z 153 Th) of quercetin aglycone was the result of a Retro-Diels-Alder (RDA) fragmentation, while the second fragment corresponded to the [M+H-H<sub>2</sub>O-2CO]<sup>+</sup> ion.

Flavan-3-ols were mainly characterized by fragments resulting from RDA fragmentation (-152 Th), oligomeric and ethyl-bridged proanthocyanidins being also characterized by the loss of monomeric entities.

Taxifolin was qualified by its RDA fragment at m/z 153 Th, and astilbin by the loss of its rhamnose entity (-146 Th). Anthocyanin detection was mostly based on the loss of the (acetyl-/coumaroyl-/caffeoyl-)glucose entity (acetyl-glucose: -204 Th loss, coumaroyl-glucose: -308 Th loss, caffeoyl-glucose: -324 Th loss, glucose: -162 Th loss), anthocyanins 3,5-diglucosides undergoing the loss of two glucose entities. (epi)catechin-(ethyl-)anthocyanins adducts were also characterized by the additional loss of the (epi)catechin entity.

**Table 1.** List of quantified compounds, MRM parameters (ion mode, precursor and product ions m/z, retention times), and calibration ranges.

		I. M. I.	m/z Precursor	m/z Quantifier <sup>a</sup>	m/z Qualifiers <sup>a</sup>	Retention	Calibration
		Ion Mode	Ion <sup>a</sup> (Th)	(Th)	(Th)	Time (min)	Range (mg·L <sup>-1</sup> )
	Gallic acid	-	168.9	125.0	79.0	1.58	0.01-90
	Protocatechuic acid	-	152.9	109.0	80.9	3.08	0.01-20
Benzoic acids and	Syringic acid	+	199.0	140.1	155.1	9.98	0.01-20
derivatives	Vanillic acid	+	169.0	125.0	93.0	7.29	0.01-20
	Ethyl gallate	-	197.0	124.0	169.0	10.65	0.01-20
	Ethyl protocatechuate	-	181.0	108.0	153.0	13.43	0.01-50
	Caffeic acid	-	179.0	135.0	79.0	7.69	0.01-50
	p-coumaric acid	-	163.0	119.0	91.0	10.71	0.01-20
	Ferulic acid	-	193.0	134.0	178.0; 149.0	12.29	0.01-20
	Caftaric acid (cis- and trans- isomers)	-	311.0	149.0	179.0; 135.0	4.30; 4.94	0.01-100
	Coutaric acid (cis- and trans- isomers)		295.1	163.1	149.1; 119.1	( 40, ( 01	as equivalents of
Hydroxycinnamic acids	Coutaric acia (cis- ana trans- isomers)	-	293.1	103.1	149.1; 119.1	6.48; 6.81	trans-caftaric acid
and derivatives	Fertaric acid		325.1	193.0 134.0; 149.2	9.25	as equivalents of	
and derivatives	reriaric acia	-	323.1	193.0	134.0, 149.2	9.23	trans-caftaric acid
	GRP (cis- and trans- isomers)		616.0	149.1	167.0	6.13; 7.52	as equivalents of
	GRF (cis- una trans- isomers)	-	010.0	149.1	107.0	0.13, 7.32	trans-caftaric acid
	Ethyl caffeate	-	207.1	135.0	179.0; 133.0	18.97	0.01-30
	Ethyl coumarate		191.0	117.1	163.0	22.35	as equivalents of
	Etnyi coumarate	<u>-</u>	191.0	117.1	103.0	22.33	Ethyl caffeate
	trans-resveratrol	+	229.0	135.0	107.1	15.43	0.01-20
	cis-resveratrol	+	229.0	135.0 + 107.1 b	135.0; 107.1	18.22	as equivalents of
	Cis-resveruiroi	'	449.U	133.0 + 107.1	133.0, 107.1	10.22	trans-resveratrol
Stilbenes	viniferin (2 isomers)	+	455.1	361.2	107.1	13.92; 15.33	as equivalents of
Julioches	vingerin (2 isomers)	•	<b>⊤</b> ЈЈ.1	501.2	107.1	13.72, 13.33	trans-resveratrol
	trans-piceid	+	391.1	229.1	135.0	12.87	0.01-20
	cis-piceid	+	391.1	229.1 + 135.0 b	229.1; 135.0	15.57	as equivalents of
	cio piccia	•	5/1.1	227.1 133.0	227.1, 133.0	10.01	trans-piceid

 Table 1. Cont.

		Ion Mode	m/z Precursor	m/z Quantifier <sup>a</sup>	m/z Qualifiers <sup>a</sup>	Retention	Calibration
		ion Mode	Ion <sup>a</sup> (Th)	(Th)	(Th)	Time (min)	Range (mg·L <sup>-1</sup> )
	Quercetin Glc <sup>c</sup>	+	465.2	303.0	84.9	15.54	0.01-30
	Myricetin Glc	+	481.2	319.1	123.1	13.52	as equivalents of
	Myriceun Gic	т	461.2	319.1	123.1	13.32	Quercetin Glc
Flavonols	Myricetin glucuronide	+	495.2	319.1	84.9	13.34	as equivalents of
Flavoliois	Myriceun giucuroniae	Т	493.2	319.1	04.9	13.34	Quercetin Glc
	Quercetin glucuronide	+	479.2	303.1	85.1; 153.1	15.21	as equivalents of
		Т	479.2	303.1	65.1, 155.1		Quercetin Glc
	Quercetin	+	303.0	153.1	229.1	21.51	0.01-10
	Catechin	+	291.1	139.0	123.1; 165.1	6.95	0.01-50
	Galloylated dimer	+	731.2	127.0	139.0; 123.0	11.29	as equivalents of Catechin
	Trimer-1 (1 isomer)	+	867.3	127.1	579.2; 289.1	3.47	as equivalents of Catechin
	Trimers-2 (several isomers)	+	867.3	127.1	579.2; 289.1	5.5-14.0	as equivalents of Catechin
	(epi)cat-ethyl-(epi)cat-1 (1 isomer) <sup>d</sup>	+	607.1	317.2	165.1	13.91	as equivalents of Catechin
Flavan-3-ols	(epi)cat-ethyl-(epi)cat-2 (2 co-eluted isomers)	+	607.1	317.2	165.1	21.37	as equivalents of Catechin
	Epicatechin	+	291.1	139.1	123.1; 165.1	10.84	0.01-50
	Dimer B2	+	579.2	127.0	139.1; 289.1	9.89	0.01-15
	Dimer B1	+	579.2	127.0	139.1; 289.1	6.74	as equivalents of Dimer B2
	Dimer B3	+	579.3	127.1	139.1; 289.2	6.21	as equivalents of Dimer B3
	Dimer B4	+	579.2	127.0	139.1; 289.1	8.20	as equivalents of Dimer B2
	A. dill. iv.		451.1	205.1	95 0. 147 1	14.10	as equivalents of
Dihydroflavonols	Astilbin	+	451.1	305.1	85.0; 147.1	14.18	Quercetin Glc
	Taxifolin	+	305.0	259.1	153.0	12.28	0.01-20

 Table 1. Cont.

		Ion Mode	m/z Precursor	m/z Quantifier <sup>a</sup>	m/z Qualifiers <sup>a</sup>	Retention	Calibration
		ion wiode	Ion <sup>a</sup> (Th)	(Th)	(Th)	Time (min)	Range (mg·L <sup>-1</sup> )
	Malvidin 3,5-diGlc	+	655.2	331.2	493.2	11.87	0.01-30
	Delphinidin 3,5-diGlc	+	627.2	303.2	465.2	9.74	as equivalents of
	Detprimain 3,5-artic	1	027.2	303.2	403.2	9.7 <b>4</b>	Malvidin 3,5-diGlc
	Cyanidin 3,5-diGlc	+	611.2	287.2	449.2	10.71	as equivalents of
	Cyaniain 5,5-aiGic	1	011.2	207.2	449.2	10.71	Malvidin 3,5-diGlc
	Petunidin 3,5-diGlc	+	641.2	317.2	479.2	11.13	as equivalents of
	1 etumum 3,5-atote	'	041.2	317.2	7/2.2	11.13	Malvidin 3,5-diGlc
	Peonidin 3,5-diGlc	+	625.2	301.2	463.2	11.70	as equivalents of
	1 contain 3,3-aiGic	'	023.2	301.2	403.2	11.70	Malvidin 3,5-diGlc
	Malvidin 3-O-Glc	+	493.2	331.2	315.0	13.80	0.01-50
	Delphinidin 3-O-Glc	+	465.2	303.1	229.0	11.78	as equivalents of
	Delphinain 3-0-01c	'	403.2	303.1	227.0	11.70	Malvidin 3-O-Glc
nthocyanins	Cyanidin 3-O-Glc	+	449.2	287.1	137.2	12.37	as equivalents of
ntnocyannis	Cyamam 5-0-010	'	777.2	207.1	137.2	12.57	Malvidin 3-O-Glc
	Petunidin 3-O-Glc	+	479.2	317.2	302.3	12.87	as equivalents of
	1 etaman 5-0-0te	'	₹/).2	317.2	302.3	12.07	Malvidin 3-O-Glc
	Peonidin 3-O-Glc	+	463.2	301.2	286.1	13.40	as equivalents of
	1 contain 5-0-dic	'	403.2	301.2	200.1	13.40	Malvidin 3-O-Glc
	Delphinidin 3-O-acetyl-Glc	+	507.2	303.0	229.0	16.04	as equivalents of
	Delphinian 3-0-acetyl-Gic	'	307.2	303.0	227.0	10.04	Malvidin 3-O-Glc
	Cyanidin 3-O-acetyl-Glc	+	491.2	287.1	137.0	18.33	as equivalents of
	Cyamam 5-0-acetyr-Ote	'	771.2	207.1	137.0	10.55	Malvidin 3-O-Glc
	Petunidin 3-O-acetyl-Glc	+	521.2	317.2	302.3	19.76	as equivalents of
	1 cianum 5-0-uceiyi-0ic	1	321.2	31/.2	302.3	17.70	Malvidin 3-O-Glc
	Peonidin 3-O-acetyl-Glc	+	505.2	301.1	286.2	20.97	as equivalents of
	1 comun 5-0-uceryr-Oic	ı	505.2	501.1	200.2	20.71	Malvidin 3-O-Glc

 Table 1. Cont.

		T Mr. 1.	m/z Precursor	m/z Quantifier <sup>a</sup>	m/z Qualifiers <sup>a</sup>	Retention	Calibration
		Ion Mode	Ion a (Th)	(Th)	(Th)	Time (min)	Range (mg·L <sup>-1</sup> )
	Malvidin 3-O-acetyl-Glc	+	535.2	331.2	315.	21.18	as equivalents of
	maiviain 5-0-aceiyi-Gic	ı	333.2	331.2	313.	21.10	Malvidin 3-O-Glc
	Delphinidin 3-O-coumaroyl-Glc	+	611.2	303.1	229.1	20.71; 22.05	as equivalents of
	(cis- and trans- isomers)	Т	011.2	303.1	229.1	20.71, 22.03	Malvidin 3-O-Glc
	Cyanidin 3-O-coumaroyl-Glc	+	595.2	287.1	157.1	21.60; 22.67	as equivalents of
	(cis- and trans- isomers)	ı	393.2	207.1	137.1	21.00, 22.07	Malvidin 3-O-Glc
	Petunidin 3-O-coumaroyl-Glc	+	625.2	317.2	302.3	21.97; 22.83	as equivalents of
	(cis- and trans- isomers)	1	023.2	317.2	302.3	21.77, 22.03	Malvidin 3-O-Glc
	Peonidin 3-O-coumaroyl-Glc	+	609.2	301.2	286.2	22.53; 23.17	as equivalents of
	(cis- and trans- isomers)	'	007.2	301.2	200.2	22.33, 23.17	Malvidin 3-O-Glc
	Malvidin 3-O-coumaroyl-Glc	+	639.2	331.2	315.2	22.60; 23.18	as equivalents of
	(cis- and trans- isomers)	•					Malvidin 3-O-Glc
Anthocyanins	Delphinidin 3-O-caffeoyl-Glc	+	627.2	303.1	569.2	20.21	as equivalents of
7 minocy annis	Deiphinian 5 0 eagreogr Gre	•	027.2	303.1	307.2	20.21	Malvidin 3-O-Glc
	Cyanidin 3-O-caffeoyl-Glc	+	611.2	287.0	_	21.27	as equivalents of
	eyamam v o cayyecyv ou		011. <b>2</b>	207.0			Malvidin 3- <i>O</i> -Glc
	Petunidin 3-O-caffeoyl-Glc	+	641.2	317.2	_	21.69	as equivalents of
	2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2						Malvidin 3- <i>O</i> -Glc
	Peonidin 3-O-caffeoyl-Glc	+	625.2	301.2	286.2	22.35	as equivalents of
							Malvidin 3- <i>O</i> -Glc
	Malvidin 3-O-caffeoyl-Glc	+	655.2	331.2	315.1	22.43	as equivalents of
							Malvidin 3- <i>O</i> -Glc
	(epi)cat-ethyl-peonidin 3-O-Glc (4 isomers)	+	779.2	327.1	489.1	16.98; 18.38;	as equivalents of
	(1)				.02.2	20.26; 20.90	Malvidin 3- <i>O</i> -Glc
	(epi)cat-ethyl-malvidin 3-O-Glc (4 isomers)	+	809.2	357.3	519.1	17.93; 19.46; 20.46	as equivalents of
-	(-1-7				- /	(2 co-eluted isomers)	Malvidin 3- <i>O</i> -Glc

 Table 1. Cont.

		Ion Mode	m/z Precursor Ion <sup>a</sup> (Th)	m/z Quantifier <sup>a</sup> (Th)	m/z Qualifiers <sup>a</sup> (Th)	Retention Time (min)	Calibration Range (mg·L <sup>-1</sup> )
	(epi)cat-ethyl-malvidin 3-O-coumaroyl-Glc (2 co-eluted isomers)	+	955.2	647.2	357.1	23.31	as equivalents of Malvidin 3-O-Glc
	Delphinidin 3-O-Glc-(epi)cat	+	753.2	591.2	303.1	8.32	as equivalents of Malvidin 3- <i>O</i> -Glc
	Cyanidin 3-O-Glc-(epi)cat	+	737.2	575.2	287.1	8.53	as equivalents of Malvidin 3- <i>O</i> -Glc
	Petunidin 3-O-Glc-(epi)cat	+	767.2	605.2	317.1	10.20	as equivalents of Malvidin 3- <i>O</i> -Glc
	Peonidin 3-O-Glc-(epi)cat	+	751.2	589.2	301.1	10.98	as equivalents of Malvidin 3-O-Glc
	Malvidin 3-O-Glc-(epi)cat (2 isomers)	+	781.2	619.2	257.1	11.05; 12.09	as equivalents of Malvidin 3-O-Glc
Anthocyanins	Malvidin 3-O-coumaroyl-Glc-(epi)cat (2 isomers)	+	927.2	619.2	373.2	15.43; 15.96	as equivalents of Malvidin 3-O-Glc
	(epi)cat-delphinidin 3-O-Glc A-F bicyclic	+	755.2	315.2	593.2	11.13	as equivalents of Malvidin 3-O-Glc
	(epi)cat-cyanidin 3-O-Glc A-F bicyclic	+	739.2	299.1	587.2	11.78	as equivalents of Malvidin 3-O-Glc
	(epi)cat-petunidin 3-O-Glc A-F bicyclic	+	769.2	329.1	617.2	12.20	as equivalents of Malvidin 3-O-Glc
	(epi)cat-peonidin 3-O-Glc A-F bicyclic	+	753.2	313.1	601.2	12.87	as equivalents of  Malvidin 3-O-Glc
	(epi)cat-malvidin 3-O-Glc A-F bicyclic	+	783.2	343.1	469.2	13.14	as equivalents of
	Caftaric-peonidin 3-O-Glc (2 co-eluted isomers)	+	773.2	611.2	151.0	11.07	Malvidin 3- <i>O</i> -Glc as equivalents of Malvidin 3- <i>O</i> -Glc

 Table 1. Cont.

		Ion Mode	m/z Precursor	m/z Quantifier <sup>a</sup>	m/z Qualifiers <sup>a</sup>	Retention	Calibration
		ion wrode	Ion <sup>a</sup> (Th)	(Th)	(Th)	Time (min)	Range (mg·L <sup>-1</sup> )
	Caftaric-malvidin 3-O-Glc (2 co-eluted isomers)	+	803.2 641.2	641.2	181.1	11.62	as equivalents of
	Cajiaric-maiviain 5-0-Gic (2 co-eiulea isomers)	Т	803.2	041.2	101.1	11.02	Malvidin 3-O-Glc
	Contania mahidin 2 O Cla (2 as alutad isamana)	1	787.2	625.2		12.44	as equivalents of
	Coutaric-malvidin 3-O-Glc (2 co-eluted isomers)	+	101.2	023.2	-	12.44	Malvidin 3-O-Glc
	Pyranodelphinidin 3-O-Glc	+	489.2	327.2		12.19	as equivalents of
	r yrunoueipniniain 3-0-0ic	Т	469.2	321.2	-	12.19	Malvidin 3-O-Glc
	Pyranocyanidin 3-O-Glc	+	473.2	311.2		13.04	as equivalents of
	1 yrunocyuniain 5-0-0ic	1	473.2	311.2	-	13.04	Malvidin 3-O-Glc
	Pyranopetunidin 3-O-Glc	+	503.2	341.1		13.69	as equivalents of
	1 yrunopeiunium 3-0-0ic	1	303.2	341.1		13.09	Malvidin 3-O-Glc
	Pyranopeonidin 3-O-Glc	+	487.2	325.1	310.2	14.67	as equivalents of
	1 yrunopeomum 5-0-01c	1	467.2	323.1	310.2	14.07	Malvidin 3-O-Glc
Anthocyanins	Pyranomalvidin 3-O-Glc (vitisin B)	+	517.2	355.1	339.2	15.55	as equivalents of
Anthocyannis	1 yrunomatviain 5-0-0tc (vitisin B)	1	317.2	333.1	339.2	13.33	Malvidin 3-O-Glc
	Carboxypyranodelphinidin 3-O-Glc	+	533.2	371.2		12.36	as equivalents of
	Сигоохуругиной гринийн 3-0-010	1	333.2	371.2	-	12.30	Malvidin 3-O-Glc
	Carboxypyranocyanidin 3-O-Glc	+	517.2	355.2	_	13.26	as equivalents of
	Carboxypyranocyanian 3-0-0ic	'	317.2	333.2	-	13.20	Malvidin 3-O-Glc
	Carboxypyranopetunidin 3-O-Glc	+	547.2	385.2	271.2	13.80	as equivalents of
	Carboxypyranoperanium 3-0-0tc		547.2	363.2	271.2	13.00	Malvidin 3-O-Glc
	Carboxypyranopeonidin 3-O-Glc	+	531.2	369.2	297.2	14.97	as equivalents of
	Carooxypyrunopeonium 5-0-0ic		331.2	507.2	<u> </u>	11.77	Malvidin 3-O-Glc
	Carboxypyranomalvidin 3-O-Glc (vitisin A)	+	561.2	399.1	383.1	15.87	as equivalents of
	Carooxypyrunomaviain 5-0-0ic (viiisii A)	₹	301.2	377.1	303.1	15.07	Malvidin 3- <i>O</i> -Glc
	Pyranodelphinidin 3-O-acetyl-Glc	+	531.2	327.2	_	12.79	as equivalents of
	1 угановегринши 5-0-исегуг-ос	1	JJ1.4	341.4	_	14.17	Malvidin 3-O-Glc

 Table 1. Cont.

		Ion Mode	m/z Precursor	m/z Quantifier <sup>a</sup>	m/z Qualifiers <sup>a</sup>	Retention	Calibration
		Ion Mode	Ion <sup>a</sup> (Th)	(Th)	(Th)	Time (min)	Range (mg·L <sup>-1</sup> )
	Pyranocyanidin 3-O-acetyl-Glc	+	515.2	311.2		13.91	as equivalents of
	r yranocyaniain 3-0-aceiyi-Gic	Т	313.2	311.2	-	13.91	Malvidin 3-O-Glc
	Pyranopetunidin 3-O-acetyl-Glc	+	545.2	341.2		14.50	as equivalents of
	r yrunopetunium 5-0-ucetyt-Gic	Т	343.2	341.2	-	14.30	Malvidin 3-O-Glc
	Pyranopeonidin 3-O-acetyl-Glc	+	529.2	325.2	309.2	15.97	as equivalents of
	1 yrunopeomum 3-0-ucetyt-Gic	1	329.2	323.2	309.2	13.97	Malvidin 3-O-Glc
	Pyranomalvidin 3-O-acetyl-Glc	+	559.2	355.2	339.2	16.67	as equivalents of
	1 yrunomatviain 3-0-acetyi-Gic	1	337.2	333.2	337.2	10.07	Malvidin 3-O-Glc
	Carboxypyranodelphinidin 3-O-acetyl-Glc	+	575.2	371.2	283.2	12.62	as equivalents of
	Carboxypyranoaeipniniain 5-0-aceiyi-0ic	ı	373.2	3/1.2	263.2	12.02	Malvidin 3-O-Glc
	Carboxypyranocyanidin 3-O-acetyl-Glc	+	559.2	355.2	527.2	13.91	as equivalents of Malvidin
	Carooxypyranocyaniain 5-0-acetyi-Gic		337.2	333.2	327.2	13.71	3-O-Glc
Anthocyanins	Carboxypyranopetunidin 3-O-acetyl-Glc	+	589.2	385.2	_	14.34	as equivalents of
Anthocyannis	Carboxypyranopeianiain 3-0-acetyt-01c	ı	389.2	363.2	-	14.34	Malvidin 3-O-Glc
	Carboxypyranopeonidin 3-O-acetyl-Glc	+	573.2	369.2	281.2	16.28	as equivalents of
	Carboxypyranopeoniain 5-0-acetyt-Gic	1	373.2	307.2	201.2	10.26	Malvidin 3-O-Glc
	Carboxypyranomalvidin 3-O-acetyl-Glc	+	603.2	399.2	383.2	16.90	as equivalents of
	Carboxypyranomarviain 5-0-acetyt-Gic	1	003.2	377.2	363.2	10.50	Malvidin 3-O-Glc
	Pyranopeonidin 3-O-coumaroyl-Glc	+	633.2	325.2	310.2	21.69	as equivalents of
	1 yrunopeonium 5-0-coumuroyi-Gic	1	033.2	323.2	310.2	21.07	Malvidin 3-O-Glc
	Pyranomalvidin 3-O-coumaroyl-Glc	+	663.2	355.2		21.78	as equivalents of
	1 yrunomurvium 3-0-coumuroyi-0ic	1	003.2	333.2	-	21.76	Malvidin 3-O-Glc
	Carboxypyranopetunidin 3-O-coumaroyl-Glc	+	693.2	385.2	370.2	19.13	as equivalents of
	Сагоохуруганорешнийн 3-0-соинагоун-онс	1	073.2	303.2	370.2	17.13	Malvidin 3-O-Glc
	Carboxypyranopeonidin 3-O-coumaroyl-Glc	+	677.2	369.2	354.1	21.37	as equivalents of
	сагоохуруганореонийн 3-0-сойтагоуг-ӨС	ı	011.2	307.2	JJ4.1	41.3/	Malvidin 3-O-Glc

 Table 1. Cont.

		Ion Mode	m/z Precursor	m/z Quantifier <sup>a</sup>	m/z Qualifiers <sup>a</sup>	Retention	Calibration
		Ion Mode	Ion <sup>a</sup> (Th)	(Th)	(Th)	Time (min)	Range (mg·L <sup>-1</sup> )
	Carboxypyranomalvidin 3-O-coumaroyl-Glc	+	707.2	399.0	383.0	21.44	as equivalents of
	Carvoxypyranomaiviain 3-O-coumaroyi-Gic	Т	707.2	399.0	363.0	21.44	Malvidin 3-O-Glc
	p-hydroxyphenylpyranopeonidin 3-O-Glc	+	579.2	417.2	402.2	23.52	as equivalents of
	p-nyaroxypnenyipyranopeoniain 3-0-Gic	Т	319.2	417.2	402.2	23.32	Malvidin 3-O-Glc
	p-hydroxyphenylpyranomalvidin 3-O-Glc	+	609.2	447.2	431.1	23.56	as equivalents of
	p-nyaroxypnenyipyranomaiviain 3-0-Gic	Т	009.2	447.2	431.1	23.30	Malvidin 3-O-Glc
	p-hydroxyphenylpyranopeonidin 3-O-acetyl-Glc	+	621.2	417.2		23.59	as equivalents of
	p-nyaroxypnenyipyranopeoniain 3-0-aceiyi-Gic	Т	621.2	417.2	-	23.39	Malvidin 3-O-Glc
	p-hydroxyphenylpyranomalvidin 3-O-acetyl-Glc	+	651.2	447.2	431.2	23.63	as equivalents of
	p-nyaroxypnenyipyranomaiviain 3-0-aceiyi-Gic	Т	031.2	447.2	431.2	23.03	Malvidin 3-O-Glc
	p-hydroxyphenylpyranopeonidin	+	725.2	417.2	_	23.67	as equivalents of
	3-O-coumaroyl-Glc	Т	123.2	417.2	-	23.07	Malvidin 3-O-Glc
nthocyanins	p-hydroxyphenylpyranomalvidin	+	755.2	447.1	431.1	23.68	as equivalents of
minocyanins	3-O-coumaroyl-Glc	Т	133.2	44 / . 1	431.1	23.08	Malvidin 3-O-Glc
	Catechylpyranopeonidin 3-O-Glc	+	595.2	433.2	418.1	23.39	as equivalents of
	Catecnytpyranopeoniain 5-0-Gic	T	393.2	433.2	410.1	23.39	Malvidin 3-O-Glc
	Catechylpyranomalvidin 3-O-Glc (pinotin A)	+	625.2	463.2	447.2	23.44	as equivalents of
	Catecnytpyranomatviain 5-O-Gic (pinotin A)	Т	623.2	403.2	447.2	23.44	Malvidin 3-O-Glc
	Catechylpyranopetunidin 3-O-acetyl-Glc	+	653.2	449.2		23.59	as equivalents of
	Catecnytpyranopetuniath 3-0-acetyt-Gic	Т	033.2	449.2	-	23.39	Malvidin 3-O-Glc
	Catachuluman an aonidin 2 O acetul Cla	1	637.2	433.2		23.47	as equivalents of
	Catechylpyranopeonidin 3-O-acetyl-Glc	+	037.2	433.2	-	23.47	Malvidin 3-O-Glc
	Catechylpyranomalvidin 3-O-acetyl-Glc	_	667.2	463.2		22 51	as equivalents of
	Catecnytpyranomatviain 3-O-acetyt-Glc	+	00/.2	403.2	-	23.51	Malvidin 3-O-Glc
	Cataoladanaga an atauridin 2 O acum Cl-	1	757.3	440.2		22.67	as equivalents of
	Catechylpyranopetunidin 3-O-coumaroyl-Glc	+	757.2	449.2	-	23.67	Malvidin 3-O-Glc

 Table 1. Cont.

		Ion Mode	m/z Precursor	m/z Quantifier <sup>a</sup>	m/z Qualifiers <sup>a</sup>	Retention	Calibration
		Ion Mode	Ion <sup>a</sup> (Th)	(Th)	(Th)	Time (min)	Range (mg·L <sup>-1</sup> )
	Catechylpyranopeonidin 3-O-coumaroyl-Glc	+	741.2	433.2		23.55	as equivalents of
	Catecnyipyranopeontain 5-0-coumaroyi-Gic	Т	/41.2	433.2	-	23.33	Malvidin 3-O-Glc
	Catechylpyranomalvidin 3-O-coumaroyl-Glc	+	771.2	463.2		23.59	as equivalents of
	Catecnyipyranomatviain 3-0-coumaroyi-Gic	Т	//1.2	403.2	-	23.39	Malvidin 3-O-Glc
	Guaiacylpyranomalvidin 3-O-Glc	+	639.2	477.2	462.2	23.60	as equivalents of
	Guaiacyipyrunomaiviain 3-0-01c	ı	039.2	4//.2	402.2	23.00	Malvidin 3-O-Glc
	Guaiacylpyranomalvidin 3-O-acetyl-Glc	+	681.2	477.2	_	23.63	as equivalents of
	Guatacytpyranomatviain 3-0-acetyt-Gic	1	001.2	7//.2	-	23.03	Malvidin 3-O-Glc
	Guaiacylpyranomalvidin 3-O-coumaroyl-Glc	+	785.2	477.1	461.3	23.71	as equivalents of
	Guatacytpyranomatviain 3-0-coumaroyi-Gic	1	763.2	7//.1	401.5	23.71	Malvidin 3-O-Glc
	Syringylpyranomalvidin 3-O-Glc	+	669.2	507.2 -	_	23.59	as equivalents of
	Syringyipyrunomawain 5-0-01c	'	007.2				Malvidin 3-O-Glc
Anthocyanins	Petunidin 3-O-acetyl-Glc-(epi)cat	+	809.2	359.1	605.2; 317.1	13.05	as equivalents of
7 minoc y annis	Teuman 5 0 accept the (epipear	•	007.2	337.1	003.2, 317.1	13.03	Malvidin 3-O-Glc
	Peonidin 3-O-acetyl-Glc-(epi)cat	+	793.2	343.1	589.2; 301.1	13.22	as equivalents of
	Teomain's o access one (epi)ear	•	175.2	313.1	307.2, 301.1	13.22	Malvidin 3-O-Glc
	Malvidin 3-O-acetyl-Glc-(epi)cat	+	823.2	373.2	331.1; 619.2	13.49	as equivalents of
	marram 5 8 accept the (epi)car	•	023.2	313.2	331.1, 017.2	13.17	Malvidin 3-O-Glc
	Pyranopeonidin 3-O-Glc-(epi)cat	+	775.2	613.2	461.2	23.19	as equivalents of
	Tyrunopeomain 5 0 Gie (epiyeur		773.2	013.2	101.2	23.19	Malvidin 3-O-Glc
	Pyranomalvidin 3-O-Glc-(epi)cat	+	805.2	491.0	643.2	23.43	as equivalents of
	Tyrunomarriam 5 0 Ole (epiyeur		000.2	1,71.0	0.13.2	23.13	Malvidin 3-O-Glc
	Pyranomalvidin 3-O-coumaroyl-Glc-(epi)cat	+	951.2	643.2	_	23.43	as equivalents of
	- y. aoa 2 0 communoyi Gio (opi)cui	ı		3 ·2·-		25.75	Malvidin 3-O-Glc
	Unknown "551" e	+	551.2	221.2	151.1	13.8	as equivalents of
	Charlette DD1	•	221.2	1.2	131.1	15.0	Malvidin 3-O-Glc

Table 1. Cont.

		T Mr. J.	m/z Precursor	m/z Quantifier <sup>a</sup>	m/z Qualifiers <sup>a</sup>	Retention	Calibration
		Ion Mode	Ion <sup>a</sup> (Th)	(Th)	(Th)	Time (min)	Range (mg·L <sup>-1</sup> )
	11.1		591.2	251.2	165.1	1.4.42	as equivalents of
	Unknown "581" <sup>e</sup>	+	581.2	251.2	165.1	14.43	Malvidin 3-O-Glc
	Unknown "1095" e		1005.0	771 0	202.1	22.25	as equivalents of
Anthocyanins		+	1095.2	771.2	302.1	22.35	Malvidin 3-O-Glc
	Unknown "1099" °	+	1000.2	027.2	- 21.18	21.10	as equivalents of
			1099.2	937.2		21.18	Malvidin 3-O-Glc
	Tyrosine	+	182.0	136.1	165.1	1.61	0.01-50
	Tyrosol <sup>f</sup>	+	121.0	77.0	93.1	5.3	0.01-50
Others (aminoacids	Hydroxytyrosol	+	155.0	91.0	81.0	3.61	0.01-20
and alcohols)	Tryptophol	+	162.0	127.2	117.0	12.36	0.01-20
	Tryptophan	-	203.0	116.0	74.0	5.98	0.01-50

<sup>&</sup>lt;sup>a</sup> span: 0.2; <sup>b</sup> The *cis*-isomers of resveratrol and piceid were quantified by integrating a sum of two transition signals, as explained in part 2.1; <sup>c</sup> Glc: glucoside; <sup>d</sup> epi(cat): catechin or epicatechin; <sup>e</sup> The reasons for the quantification of the unknown compounds as equivalents of malvidin 3-*O*-glucoside are presented in part 2.2; <sup>f</sup> The precursor ion *m/z* for tyrosol does not refer to the [M+H]<sup>+</sup> ion, but to the [M-H<sub>2</sub>O+H]<sup>+</sup> ion due to in-source fragmentation, even under mild ionization conditions [45].

Tyrosol was characterized by two fragments: m/z 77 Th and m/z 93 Th corresponding to the aromatic ring and the phenol group, respectively. Tryptophan's detection was based on the fragment at m/z 116 Th representing the indole ion of the molecule obtained by the loss of the side-chain, and the fragment at m/z 74 Th representing the [CHNH<sub>2</sub>-COO]<sup>-</sup> ion. Tyrosine was characterized by the formation of an immonium ion and by the loss of a -NH<sub>3</sub> group.

Polyphenols are a wide family of natural compounds, therefore the major part of them are not commercially available, and their deuterated forms even less so. For that reason, an absolute quantification was not possible, and these compounds were quantified as equivalents of reference standards from the same chemical class, as specified in Table 1. Such a relative quantification method is a suitable approach since the current UHPLC-QqQ-MS method is intended to large-scale comparative studies. Several particular cases were pointed out. The cis-resveratrol fragmentation pathway is quantitatively different from the one of the trans-isomer. Quantifying cis-resveratrol as equivalents of trans-resveratrol, based on the major transition (or even a minor one), was thus not appropriate, by comparison to UV-based quantification [47]. It appeared that quantifying cis-resveratrol as equivalents of trans-resveratrol by integrating the sum of both selected transition traces instead of only one of them matched to the UV-based quantification. This method was verified on several wine samples spiked with different amounts of resveratrol (data not shown). This statement was also applied by extension to the quantification of cis-piceid (the glucoside form of resveratrol) as equivalents of trans-piceid. Ethyl coumarate was not commercially available and was quantified as equivalents of ethyl caffeate. However, this quantification was not based on the main transition of ethyl caffeate (Electrospray Ionization (ESI) ESI + 207 > 135 Th) but on a minor one (ESI + 207 > 133 Th; -74 Th loss) so as to match to the main fragmentation of ethyl coumarate (ESI + 191 > 117 Th; -73 Th loss). In the same way, coutaric acid (ESI - 295 > 163 Th) and fertaric acid (ESI - 325 > 193 Th) quantifications were based on a minor, still intense, fragmentation transition of caftaric acid (ESI -311 > 179 Th) corresponding to the same -132 Th loss. GRP (main fragmentation: ESI - 616 > 149 Th) was quantified by resorting to the major fragmentation of caftaric acid (ESI -311 > 149 Th) in order to obtain the same fragment. Astilbin, also known as dihydroquercetin 3-O-rhamnoside, is not commercially available, and was thus quantified as equivalents of its structurally closest compound, quercetin glucoside, which also loses its sugar moiety.

## 2.2. Unknown Compounds

Four compounds detected by UHPLC-IT-MS were not identified. These compounds were named according to their m/z ratio in positive mode: "unknown581" and "unknown551", both consisting in double peaks, "unknown1099" which was detected as doubly charged [550]<sup>2+</sup>, and "unknown1095". They displayed an absorbance in the anthocyanin wavelength range (around 520 nm) and were characterized by the loss of one glucose entity (-162) for "unknown581" and "unknown551" and two successive glucose entity losses for "unknown1099" and "unknown1095". The m/z difference between "unknown 551" and "unknown 581" was also noticeable, since it corresponded to the mass difference between peonidin and malvidin compounds. Although their structure was not elucidated, these unknown compounds were added to the UHPLC-QqQ-MS method because of their potential relevancy for metabolomics' approaches. They were quantified as equivalents of malvidin 3-*O*-glucoside, as they displayed UV spectra and

glucose losses characteristic of anthocyanin compounds. Fragmentation spectra are available as supplementary data.

# 2.3. Optimisation of UHPLC-QqQ-MS Conditions

UHPLC is a technique that allows an enhanced chromatographic resolution compared to classical HPLC, due to an increase in column pressure (up to 1000 bars). The optimal working flow rate is then much higher with UHPLC than with HPLC. This consequence motivated the choice of a 1 mm-diameter column versus a 2.1 mm-diameter column for the development of the quantification method described in this paper. The mobile phase flow rate could thus be decreased from 0.5-0.6 mL·min<sup>-1</sup> to 0.17 mL·min<sup>-1</sup> in order to suit better to the use of an ESI-MS interface, without degrading the peak resolution, while reducing solvent consumption. The percentage of formic acid is a key parameter for polyphenol analysis, since it has to be high enough for maintaining the anthocyanins in their flavylium cationic form (pH  $\leq$  2) [48], but not too high so as to avoid signal suppression in the source of the mass spectrometer, especially for acids. Three percentage values (1%, 2% and 5%) were tested for the analysis of several hydroxycinnamic acids. A MS signal intensity loss of 25% and 50% was observed when increasing the percentage of formic acid from 1% to 2% and 5% respectively. The optimal working percentage of formic acid in the mobile phase was thus set at 1%, what corresponded to a pH of 2.1 in solvent A. Methanol as solvent B was preferred to acetonitrile despite its higher viscosity, because an enhancement of the MS signal intensity was observed when using methanol, especially for acids in negative ionization mode [49].


Polyphenolic compounds may be degraded by heat [50], however, increasing the column temperature to 40 °C was necessary to achieve optimal chromatographic resolution by reducing methanol viscosity and thus increasing the mobile phase flow rate.

Once the MRM transitions parameters optimised as described in Table 1, the chromatographic separation was adapted so as to obtain a gradient elution program as short as possible without losing sensitivity. Sensitivity was maintained by adapting dwell times for each compound MS transition so as to obtain between 10 and 15 points per peak. Dwell times ranged from 0.005 s, *i.e.*, the higher scanning speed limitation of the QqQ mass spectrometer used in MRM mode, to 0.705 s, with a mean and median dwell times of 0.077 s and 0.005 s respectively. The higher the number of overlapped MRM functions, the lower the dwell times. Most of the dwell times were very low, which showed the high-throughput character of the method exploiting the full capacities of the mass spectrometer analyser. The wide gap between the mean and the median values illustrated the high heterogeneity of the analysed compounds. This resulted in different elution behaviours, hence a variation in the number of overlapped MRM functions during the chromatographic separation. Some time sequences were characterized by a great number of coeluted compounds, whereas other time sequences were relatively poorer in terms of analysed compounds.

## 2.4. Advantages of the Present Method

One hundred and fifty four (154) compounds were quantified by UHPLC-QqQ-MS within a chromatographic separation of 30 min including column re-equilibration. The wine samples were analysed without prior sample extraction, clean-up or concentration steps.

Such a short analysis time could be obtained thanks to the selectivity of the method: the MRM mode focuses on transitions that are specific to the quantified compounds. Compounds which show the same precursor and product ions, *i.e.*, delphinidin 3-*O*-glucoside and quercetin glucoside or anthocyanins 3-*O*-caffeoyl-Glc and anthocyanins 3,5-diGlc, remained separated by their retention times. An example of a chromatogram is shown in Figure 1.



**Figure 1.** Superimposition of MRM traces from the UHPLC-QqQ-MS analysis of a rosé wine, (a) full display; (b) six-time zoomed display.

This example highlights the high heterogeneity of the wine sample with some major compounds and a large number of minor compounds. The large dynamic range of the QqQ mass spectrometer allowed the measurement of such concentration differences. The six-time zoomed chromatogram displays a multitude of overlapped MRM peaks, which illustrates the complexity of the sample and the ability of the MRM method to specifically distinguish co-eluted compounds. Improved selectivity is an important advantage compared to HPLC-DAD methods, for which compounds quantification is based on peak areas at selected wavelengths, requiring almost full separation of all compounds, with relatively long gradient elution programs. Nevertheless, the enhancement of selectivity brought by the QqQ-MRM detection did not allow to get rid of all other co-elution problems, in particular matrix effects and ionization competition between co-eluted molecules. For that reason, it was necessary to maintain a minimum chromatographic separation of 30 min while pushing the MRM scan rate of the mass spectrometer to its maximum.

The sensitivity of the present method was largely enhanced by the use of the MRM mode. Its efficiency—sensitivity enhancement and time saving—was increased by the possibility of running both positive and negative ion modes in a single run, contrary to the UHPLC-Q-TOF-MS method applied to red wines by Vaclavik *et al.* [23]. The sensitivity of this method was characterized by very low limits of

quantification (LOQ). The LOQ were determined by measuring the noise for each compound, and by converting 10 times these values into concentrations by means of the calibration curves. The concentrations of the different compounds in the 12 analyzed rosé wines with their calculated LOQ are available as supplementary data. The LOQ were different from a wine to another, as the matrix effects were specific to each wine.

Sensitivity combined to selectivity resulted in a major advance in wine polyphenolic composition analysis, since it allowed the quantification of 152 phenolic compounds (and two other non-phenolic wine compounds *i.e.*, tryptophan and tryptophol) *versus* 48 phenolic compounds with the method described by Gomez-Alonso *et al.* [13] and 57 phenolic compounds with the method described by Vrhovsek *et al.* [39]. This is all the more relevant that polyphenolic compounds are present in lesser quantities in rosé wines than in red wines, and thus more difficult to detect and quantify.


Nevertheless, a particular attention should be paid to *trans*-resveratrol, which is partially co-eluted with in-source fragmented *cis*-piceid (a glycosylated form of resveratrol which easily loses its glucose moiety).

## 2.5. Chemometric Analysis—An Illustrative Application of this Method

The present UHPLC-QqQ-MS polyphenols quantification method is particularly suitable for high-throughput wine analyses, especially when it is coupled with chemometrics. Approached in this way, large-scale food metabolomics studies can be performed in order to discriminate between wines samples or connecting their polyphenolic composition to other characteristics.

To illustrate this potential, we present in Figure 2 the Principal Component Analysis (PCA) of the data table originating from the analysis in triplicate of the 12 roses wines previously selected. Mass data were imported and treated in the Matlab environment (The Mathworks, Inc., Natick, MA, USA). PCA was performed on the variance-covariance matrix, without previous column standardization of the data table.

The plot of the PCA wines scores (Figure 2A) illustrates the variability between the polyphenolic composition of the analysed samples. In comparison, the repeatability of the method, evaluated by confidence ellipses set at 95% around the mean results is quite satisfactory. Figure 2B,C represent the loadings associated to the first and the second principal component, respectively. On these plots, the abscissa axis represents the molecular weight of each analysed compound while the ordinate axis represents its contribution to the spread out of samples on the corresponding component. Thus, for the first axis, the variability between samples of opposite sign mainly address to caftaric acid (higher concentrations measured for samples 1, 3 and 8), followed by coutaric acid, GRP and dimer B1. The second axis is more complex to interpret and is mainly characterized by the contributions of GRP (higher concentration for samples presenting a positive coordinate on the second component) and by the contribution of tyrosol, gallic acid, catechin, malvidin 3-O-glucoside, dimer B1 and malvidin 3-O-acetylglucoside (higher concentration for samples presenting a negative coordinate).



**Figure 2.** (A) PCA scores plot of the wines 1 to 12; (B) Loadings associated to the first common component A1; (C) Loadings associated to the second common component A2.

# 3. Experimental Section

## 3.1. Chemicals

Formic acid and HPLC grade methanol were purchased from VWR Prolabo (Fontenay-sous-Bois, France). Deionized water was obtained from a Milli-Q purification system (Millipore, Molsheim, France). Standards of gallic acid, vanillic acid, caffeic acid and ferulic acid were purchased from Fluka (Buchs, Switzerland). Standards of protocatechuic acid, syringic acid, p-coumaric acid, ethyl protocatechuate, tyrosol, tryptophol, tryptophan, trans-resveratrol, quercetin, catechin, epicatechin and taxifolin were purchased from Sigma-Aldrich (St. Louis, MO, USA). Standards of ethyl gallate, ethyl caffeate, tyrosine, hydroxytyrosol, malvidin 3-O-glucoside chloride and malvidin 3,5-diglucoside chloride were purchased from Extrasynthese (Geney, France). Standards of *trans*-caftaric acid, quercetin glucoside and dimer B2 were purchased from Phytolab (Vestenbergsgreuth, Germany). Standard of *trans*-piceid was purchased from Selleckchem (Houston, TX, USA). Oligomeric ProCyanidins (OPC) powder used for method development was supplied by Berkem (Gardonne, France).

#### 3.2. Samples

Twelve rosé wines, respectively numbered from 1 to 12, were selected by the Centre du Rosé (Vidauban, France) for their geographic origin diversity and their wide colour range: Mateus, Portugal (sparkling); Palatinat, Germany; Slovakia; Vaud, Switzerland; Latium, Italy; La Rioja, Spain; Provence, France (two wines); Bergerac, France; Languedoc, France; Loire, France; Champagne, France (sparkling).

# 3.3. Standard and Sample Preparation

Standard solutions were prepared in methanol/H2O 20/80 (v/v). Their calibration ranges are summarized in Table 1. Rosé wines were filtered through 0.45  $\mu$ m filters (Hydrophobic Durapore, Millipore) under vacuum in order to remove solid particles. The filtration flasks were then stoppered and the vacuum was maintained for 5 min in order to degas the wines. 1 mL of each wine was filtered through a 0.2  $\mu$ m regenerated cellulose membrane syringe filter (Phenex, Phenomenex, Le Pecq, France), and the samples were injected in triplicates for UHPLC-QqQ-MS analysis.

## 3.4. UHPLC-DAD-IT-MS<sup>n</sup> Conditions

Preliminary experiments for compound identification were carried out on an Acquity UHPLC system equipped with a diode array detector (Waters, Saint-Quentin-en-Yvelines, France) and hyphenated to an AmaZon X ion-trap mass spectrometer (Bruker Daltonique, Wissembourg, France) with electrospray ionization (ESI) operating in switching positive and negative modes. UHPLC-DAD-ESI-IT-MSn analyses were performed with a Waters ACQUITY chromatography system coupled to a Bruker Amazon X ion-trap mass spectrometer (Acquity UPLC, Waters—AmazonX IT, Bruker Daltonique). The UHPLC system included a binary pump, a cooled autosampler maintained at 7 °C and equipped with a 5-µL sample loop, a 100-μL syringe, a 30-μL needle and a DAD. The Compass HyStar software suite was used to control the instruments and to acquire and process the data. The column used for chromatographic separation was a reversed-phase Acquity HSS T3 1.8 µm 1.0 × 150 mm, (Waters) protected by a 0.2  $\mu$ m in-line filter and maintained at 40 °C. The mobile phase consisted of 1% ( $\nu/\nu$ ) formic acid in Milli-O water (solvent A) and 1% (v/v) formic acid in methanol (solvent B). Elution was performed using a gradient system at a flow rate of 0.08 mL·min<sup>-1</sup>: 2%–30% B in 15 min; 30%–50% B in 15 min, 50%-60% B in 5 min, 60%-2% B in 6min and reequilibration (2% B for 5 min). Injection volume was 0.5 μL. Detection was operated at wavelength ranging from 250 to 600 nm. ESI conditions were as follows: HV capillary voltage, 2.5 kV; drying temperature, 200 °C; Nitrogen nebulizer gas, 10 psi; Nitrogen drying gas, 5 L·min<sup>-1</sup>; capillary exit voltage, 140 V; and syringe draw rate, 10 μL·min<sup>-1</sup>. The [M-H]<sup>-</sup>, [M+H]<sup>+</sup> and [M+2H]<sup>2+</sup> ions were selected for CID fragmentation (with helium as collision gas) to produce the MS/MS spectra.

These qualitative analyses were performed on the 12 selected rosé wines in order to establish a list of quantifiable polyphenols as extensive as possible, and potentially to detect new compounds. The MRM transitions were subsequently adapted to the UHPLC-QqQ-MS method.

## 3.5. UHPLC-QqQ-MS Conditions

UHPLC-QqQ-MS analyses were carried out using an Acquity UPLC system (Waters) hyphenated to a triple quadrupole mass spectrometer (Waters) with electrospray ionization (ESI) operating in switching positive and negative modes. The UHPLC system included a binary pump, a cooled autosampler maintained at 7 °C and equipped with a 5- $\mu$ L sample loop, a 100- $\mu$ L syringe and a 30- $\mu$ L needle, and a DAD. MassLynx software was used to control the instruments and to acquire the data, and TargetLynx software was used to process the data. The column used for chromatographic separation was a reversed-phase Acquity HSS T3 1.8  $\mu$ m 1.0 × 100 mm, (Waters) protected by a 0.2  $\mu$ m in-line filter and maintained at 40 °C. The mobile phase consisted of 1% ( $\nu$ / $\nu$ ) formic acid in Milli-Q water (solvent A) and 1% ( $\nu$ / $\nu$ ) formic acid in methanol (solvent B). The source and desolvation temperatures were respectively set at 120 °C and 450 °C. Nitrogen was used as desolvation (500 L·h<sup>-1</sup>) and cone (50 L·h<sup>-1</sup>) gas. Argon was used as collision gas at a flow rate of 0.16 mL·min<sup>-1</sup>. Capillary voltage was set at 3.5 kV in positive mode and 2.8 kV in negative mode. The MRM transitions are described in Table 1.

Samples were injected into the column by using the Partial Loop with Needle Overfill injection mode, with an injection volume of 2  $\mu$ L. The following gradient elution program, at a flow rate of 0.17 mL·min<sup>-1</sup>, was performed: isocratic 1% B from 0.0 to 2.0 min, linear 1%–5% B from 2.0 to 2.1 min, linear 5%–10% B from 2.1 to 8.0 min, linear 10%–28% B from 8.0 to 12.0 min, isocratic 28% B from 12.0 to 18.0 min, linear 28%–45% B from 18.0 to 22.0 min, linear 45%–99% B from 22.0 to 23.5 min, isocratic 99% B from 23.5 to 26.5 min. At the end of this sequence, the column was brought back to initial conditions with linear 99%–1% B from 26.5 to 27.0 min, then re-equilibrated with isocratic 1% B from 27.0 to 30.0 min. Two detectors were used: DAD scanning from 210 to 600 nm (resolution 1.2 nm) and the QqQ detector.

#### 4. Conclusions

We have developed a high-throughput UHPLC-QqQ-MS analysis method for the quantification of 152 phenolic compounds, and two non-phenolic compounds in rosé wines in 30 min, without sample pre-treatment apart from filtering. Our method combines selectivity and sensitivity, both enhanced by the use of UHPLC and a QqQ-MRM acquisition mode alternating positive and negative ionization modes in a single run. The quantitative analysis of 12 rosé wines from different regions of the world showed the suitability of this method for large-scale wine analyses. Such a study is on its way and will aim to establish a typology of worldwide rosé wines by correlating chemometrically the phenolic composition to the colour, the geographic origin and the grape varieties. This quantification method can also be considered for metabolomics' approaches as well as for fraud detection (*i.e.*, illegal wine blending).

# **Supplementary Materials**

Supplementary materials can be accessed at: http://www.mdpi.com/1420-3049/20/05/7890/s1.

# Acknowledgments

The authors thank Jean-Marc Souquet and Jean-Paul Mazauric for valuable discussions and for providing relevant samples. Marine Lambert also acknowledges GIS-IBiSA (Scientific Interest Group

for Biology, Healthcare and Agronomy Infrastructures) and the Languedoc-Roussillon region for financial support.

#### **Author Contributions**

M.L., E.M. and A.V. carried out the experimental part; G. Mazerolles carried out the chemometrics analysis; J.W. helped with preliminary experiments and discussions; G. Masson, V.C. and N.S. proposed the subject; V.C. and N.S. designed and supervised the study; M.L. wrote the manuscript. All the authors read and approved the final manuscript.

#### **Conflicts of Interest**

The authors declare no conflict of interest.

#### References

- 1. Pérez-Magariño, S.; Ortega-Heras, M.; José, M.L.G.; Boger, Z. Comparative study of artificial neural network and multivariate methods to classify Spanish DO rose wines. *Talanta* **2004**, *62*, 983–990.
- 2. Puértolas, E.; Saldaña, G.; Álvarez, I.; Raso, J. Experimental design approach for the evaluation of anthocyanin content of rosé wines obtained by pulsed electric fields. Influence of temperature and time of maceration. *Food Chem.* **2011**, *126*, 1482–1487.
- 3. Wirth, J.; Caillé, S.; Souquet, J.M.; Samson, A.; Dieval, J.B.; Vidal, S.; Fulcrand, H.; Cheynier, V. Impact of post-bottling oxygen exposure on the sensory characteristics and phenolic composition of Grenache rosé wines. *Food Chem.* **2012**, *132*, 1861–1871.
- 4. Bravo, M.N.; Feliciano, R.; Silva, S.; Coelho, A.V.; Boas, L.V.; Bronze, M.R. Analysis of trans-Resveratrol analysis: Comparison of methods and contents in Muscatel fortified wines from Setúbal region in Portugal. *J. Food Compos. Anal.* **2008**, *21*, 634–643.
- 5. Buiarelli, F.; Coccioli, F.; Jasionowska, R.; Merolle, M.; Terracciano, A. Analysis of some stilbenes in Italian wines by liquid chromatography/tandem mass spectrometry. *Rapid Commun. Mass Spectrom.* **2007**, *21*, 2955–2964.
- 6. Buiarelli, F.; Coccioli, F.; Merolle, M.; Jasionowska, R.; Terracciano, A. Identification of hydroxycinnamic acid-tartaric acid esters in wine by HPLC-tandem mass spectrometry. *Food Chem.* **2010**, *123*, 827–833.
- 7. De Villiers, A.; Alberts, P.; Tredoux, A.G.J.; Nieuwoudt, H.H. Analytical techniques for wine analysis: An African perspective; a review. *Anal. Chim. Acta* **2012**, *730*, 2–23.
- 8. Dugo, P.; Cacciola, F.; Donato, P.; Airado-Rodriguez, D.; Herrero, M.; Mondello, L. Comprehensive two-dimensional liquid chromatography to quantify polyphenols in red wines. *J. Chromatogr. A* **2009**, *1216*, 7483–7487.
- 9. Fulcrand, H.; Benabdeljalil, C.; Rigaud, J.; Cheynier, V.; Moutounet, M. A new class of wine pigments generated by reaction between pyruvic acid and grape anthocyanins. *Phytochemistry* **1998**, 47, 1401–1407.
- 10. Fulcrand, H.; dos Santos, P.-J.C.; Sarni-Manchado, P.; Cheynier, V.; Favre-Bonvin, J. Structure of new anthocyanin-derived wine pigments. *J. Chem. Soc. Perkin Trans.* 1 1996, 1, 735–739.

11. Fulcrand, H.; Mane, C.; Preys, S.; Mazerolles, G.; Bouchut, C.; Mazauric, J.P.; Souquet, J.M.; Meudec, E.; Li, Y.; Cole, R.B.; *et al.* Direct mass spectrometry approaches to characterize polyphenol composition of complex samples. *Phytochemistry* **2008**, *69*, 3131–3138.

- 12. Garrido, J.; Borges, F. Wine and grape polyphenols—A chemical perspective. *Food Res. Int.* **2013**, *54*, 1844–1858.
- 13. Gómez-Alonso, S.; García-Romero, E.; Hermosín-Gutiérrez, I. HPLC analysis of diverse grape and wine phenolics using direct injection and multidetection by DAD and fluorescence. *J. Food Compos. Anal.* **2007**, *20*, 618–626.
- 14. Gruz, J.; Novák, O.; Strnad, M. Rapid analysis of phenolic acids in beverages by UPLC-MS/MS. *Food Chem.* **2008**, *111*, 789–794.
- 15. Hayasaka, Y.; Asenstorfer, R.E. Screening for potential pigments derived from anthocyanins in red wine using nanoelectrospray tandem mass spectrometry. *J. Agric. Food Chem.* **2002**, *50*, 756–761.
- 16. Jaitz, L.; Siegl, K.; Eder, R.; Rak, G.; Abranko, L.; Koellensperger, G.; Hann, S. LC-MS/MS analysis of phenols for classification of red wine according to geographic origin, grape variety and vintage. *Food Chem.* **2010**, *122*, 366–372.
- 17. Kong, Q.J.; Ren, X.Y.; Hu, N.; Sun, C.R.; Pan, Y.J. Identification of isomers of resveratrol dimer and their analogues from wine grapes by HPLC/MS(n) and HPLC/DAD-UV. *Food Chem.* **2011**, *127*, 727–734.
- 18. Landrault, N.; Larronde, F.; Delaunay, J.-C.; Castagnino, C.; Vercauteren, J.; Merillon, J.-M.; Gasc, F.; Cros, G.; Teissedre, P.-L. Levels of stilbene oligomers and astilbin in French varietal wines and in grapes during noble rot development. *J. Agric. Food Chem.* **2002**, *50*, 2046–2052.
- 19. Mazerolles, G.; Preys, S.; Bouchut, C.; Meudec, E.; Fulcrand, H.; Souquet, J.M.; Cheynier, V. Combination of several mass spectrometry ionization modes: A multiblock analysis for a rapid characterization of the red wine polyphenolic composition. *Anal. Chim. Acta* **2010**, *678*, 195–202.
- 20. Nave, F.; Teixeira, N.; Mateus, N.; de Freitas, V. The fate of flavanol-anthocyanin adducts in wines: Study of their putative reaction patterns in the presence of acetaldehyde. *Food Chem.* **2010**, *121*, 1129–1138.
- 21. Spacil, Z.; Novakova, L.; Solich, P. Analysis of phenolic compounds by high performance liquid chromatography and ultra performance liquid chromatography. *Talanta* **2008**, *76*, 189–199.
- 22. Tian, Q.; Giusti, M.; Stoner, G.; Schwartz, S. Screening for anthocyanins using high-performance liquid chromatography coupled to electrospray ionization tandem mass spectrometry with precursor-ion analysis, product-ion analysis, common-neutral-loss analysis, and selected reaction monitoring. *J. Chromatogr. A* **2005**, *1091*, 72–82.
- 23. Vaclavik, L.; Lacina, O.; Hajslova, J.; Zweigenbaum, J. The use of high performance liquid chromatography-quadrupole time-of-flight mass spectrometry coupled to advanced data mining and chemometric tools for discrimination and classification of red wines according to their variety. *Anal. Chim. Acta* **2011**, *685*, 45–51.
- 24. Baderschneider, B.; Winterhalter, P. Isolation and characterization of novel benzoates, cinnamates, flavonoids, and lignans from Riesling wine and screening for antioxidant activity. *J. Agric. Food Chem.* **2001**, *49*, 2788–2798.

25. Chiva-Blanch, G.; Urpi-Sarda, M.; Rotches-Ribalta, M.; Zamora-Ros, R.; Llorach, R.; Lamuela-Raventos, R.M.; Estruch, R.; Andres-Lacueva, C. Determination of resveratrol and piceid in beer matrices by solid-phase extraction and liquid chromatography—tandem mass spectrometry. *J. Chromatogr. A* **2011**, *1218*, 698–705.

- 26. Dehkharghanian, M.; Adenier, H.; Vijayalakshmi, M.A. Study of flavonoids in aqueous spinach extract using positive electrospray ionisation tandem quadrupole mass spectrometry. *Food Chem.* **2010**, *121*, 863–870.
- 27. Tarascou, I.; Mazauric, J.-P.; Meudec, E.; Souquet, J.-M.; Cunningham, D.; Nojeim, S.; Cheynier, V.; Fulcrand, H. Characterisation of genuine and derived cranberry proanthocyanidins by LC-ESI-MS. *Food Chem.* **2011**, *128*, 802–810.
- 28. Guillarme, D.; Casetta, C.; Bicchi, C.; Veuthey, J.L. High throughput qualitative analysis of polyphenols in tea samples by ultra-high pressure liquid chromatography coupled to UV and mass spectrometry detectors. *J. Chromatogr. A* **2010**, *1217*, 6882–6890.
- 29. Wang, D.; Lu, J.; Miao, A.; Xie, Z.; Yang, D. HPLC-DAD-ESI-MS/MS analysis of polyphenols and purine alkaloids in leaves of 22 tea cultivars in China. *J. Food Compos. Anal.* **2008**, *21*, 361–369.
- 30. Ortega, N.; Romero, M.-P.; Macià, A.; Reguant, J.; Anglès, N.; Morelló, J.-R.; Motilva, M.-J. Comparative study of UPLC-MS/MS and HPLC-MS/MS to determine procyanidins and alkaloids in cocoa samples. *J. Food Compos. Anal.* **2010**, *23*, 298–305.
- 31. Ortega, N.; Romero, M.-P.; Macià, A.; Reguant, J.; Anglès, N.; Morelló, J.-R.; Motilva, M.-J. Obtention and characterization of phenolic extracts from different cocoa sources. *J. Agric. Food Chem.* **2008**, *56*, 9621–9627.
- 32. Di Maio, I.; Esposto, S.; Taticchi, A.; Selvaggini, R.; Veneziani, G.; Urbani, S.; Servili, M. HPLC-ESI-MS investigation of tyrosol and hydroxytyrosol oxidation products in virgin olive oil. *Food Chem.* **2011**, *125*, 21–28.
- 33. Marti, M.P.; Pantaleon, A.; Rozek, A.; Soler, A.; Valls, J.; Macia, A.; Romero, M.P.; Motilva, M.J. Rapid analysis of procyanidins and anthocyanins in plasma by microelution SPE and ultra-HPLC. *J. Sep. Sci.* **2010**, *33*, 2841–2853.
- 34. Gonthier, M.-P.; Rios, L.Y.; Verny, M.-A.; Rémésy, C.; Scalbert, A. Novel liquid chromatography-electrospray ionization mass spectrometry method for the quantification in human urine of microbial aromatic acid metabolites derived from dietary polyphenols. *J. Chromatogr. B* **2003**, 789, 247–255.
- 35. Valls, J.; Millan, S.; Marti, M.P.; Borras, E.; Arola, L. Advanced separation methods of food anthocyanins, isoflavones and flavanols. *J. Chromatogr. A* **2009**, *1216*, 7143–7172.
- 36. Fabre, N.; Rustan, I.; de Hoffmann, E.; Quetin-Leclercq, J. Determination of flavone, flavonol, and flavanone aglycones by negative ion liquid chromatography electrospray ion trap mass spectrometry. *J. Am. Soc. Mass Spectrom.* **2001**, *12*, 707–715.
- 37. Arapitsas, P.; Perenzoni, D.; Nicolini, G.; Mattivi, F. Study of sangiovese wines pigment profile by UHPLC-MS/MS. *J. Agric. Food Chem.* **2012**, *60*, 10461–10471.
- 38. Ehrhardt, C.; Arapitsas, P.; Stefanini, M.; Flick, G.; Mattivi, F. Analysis of the phenolic composition of fungus-resistant grape varieties cultivated in Italy and Germany using UHPLC-MS/MS. *J. Mass Sprectrom.* **2014**, *49*, 860–869.

39. Vrhovsek, U.; Masuero, D.; Gasperotti, M.; Franceschi, P.; Caputi, L.; Viola, R.; Mattivi, F.A Versatile targeted metabolomics method for the rapid quantification of multiple classes of phenolics in fruits and beverages. *J. Agric. Food Chem.* **2012**, *60*, 8831–8840.

- 40. Rentzsch, M.; Schwarz, M.; Winterhalter, P. Pyranoanthocyanins—An overview on structures, occurrence, and pathways of formation. *Trends Food Sci. Tech.* **2007**, *18*, 526–534.
- 41. Schwarz, M.; Hofmann, G.; Winterhalter, P. Investigations on anthocyanins in wines from *Vitis vinifera* cv. Pinotage: Factors influencing the formation of Pinotin A and its correlation with wine age. *J. Agric. Food Chem.* **2004**, *52*, 498–504.
- 42. Schwarz, M.; Wabnitz, T.C.; Winterhalter, P. Pathway leading to the formation of anthocyanin-vinylphenol adducts and related pigments in red wines. *J. Agric. Food Chem.* **2003**, *51*, 3682–3687.
- 43. Nave, F.; Teixeira, N.; Mateus, N.; de Freitas, V. Hemisynthesis and structural characterization of flavanol-(4,8)-vitisins by mass spectrometry. *Rapid Commun. Mass Spectrom.* **2010**, *24*, 1964–1970.
- 44. Abad-Garcia, B.; Berrueta, L.A.; Garmon-Lobato, S.; Gallo, B.; Vicente, F. A general analytical strategy for the characterization of phenolic compounds in fruit juices by high-performance liquid chromatography with diode array detection coupled to electrospray ionization and triple quadrupole mass spectrometry. *J. Chromatogr. A* **2009**, *1216*, 5398–5415.
- 45. Gregus, P.; Vlckova, H.; Buchta, V.; Kestranek, J.; Krivcikova, L.; Novakova, L. Ultra high performance liquid chromatography tandem mass spectrometry analysis of quorum-sensing molecules of *Candida albicans*. *J. Pharm. Biomed. Anal.* **2010**, *53*, 674–681.
- 46. Pezet, R.; Perret, C.; Jean-Denis, J.B.; Tabacchi, R.; Gindro, K.; Viret, O. *delta*-Viniferin, a resveratrol dehydrodimer: One of the major stilbenes synthesized by stressed grapevine leaves. *J. Agric. Food Chem.* **2003**, *51*, 5488–5492.
- 47. Camont, L.; Cottart, C.H.; Rhayem, Y.; Nivet-Antoine, V.; Djelidi, R.; Collin, F.; Beaudeux, J.L.; Bonnefont-Rousselot, D. Simple spectrophotometric assessment of the *trans-/cis*-resveratrol ratio in aqueous solutions. *Anal. Chim. Acta* **2009**, *634*, 121–128.
- 48. De Villiers, A.; Cabooter, D.; Lynen, F.; Desmet, G.; Sandra, P. High performance liquid chromatography analysis of wine anthocyanins revisited: Effect of particle size and temperature. *J. Chromatogr. A* **2009**, *1216*, 3270–3279.
- 49. Henriksen, T.; Juhler, R.K.; Svensmark, B.; Cech, N.B. The relative influences of acidity and polarity on responsiveness of small organic molecules to analysis with negative ion electrospray ionization mass spectrometry (ESI-MS). *J. Am. Soc. Mass Spectrom.* **2005**, *16*, 446–455.
- 50. Khanal, R.C.; Howard, L.R.; Prior, R.L. Effect of heating on the stability of grape and blueberry pomace procyanidins and total anthocyanins. *Food Res. Int.* **2010**, *43*, 1464–1469.

Sample Availability: Samples of the twelve rosé wines are not available.

© 2015 by the authors; licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution license (http://creativecommons.org/licenses/by/4.0/).