


# Supplemental Materials: A Melting Curve-Based Multiplex RT-qPCR Assay for Simultaneous Detection of Four Human Coronaviruses


Zhenzhou Wan, Ya'nan Zhang, Zhixiang He, Jia Liu, Ke Lan, Yihong Hu and Chiyu Zhang

## Supplemental Results: Optimization of the Primer Concentration

In the conventional PCR, primer concentration was 0.2  $\mu$ M. When 0.2  $\mu$ M each set of primers were used to amplify a mixed template containing 30,000 copies of each HCoV RNA, only low or no melting peaks were observed for HCoV-229E, -NL63, and -HKU1, suggesting a competitive amplification of HCoV-OC43 with other three HCoVs. To obtain well performance of the melting curve-based multiplex RT-qPCR, we optimized the primer inputs for HCoV-OC43 and other three HCoVs by orthogonal experiments with primer concentration from 0.2 to 0.6  $\mu$ M. When GC-HCoV-OC43-F/R was remained at lower concentration of 0.2  $\mu$ M, the melting peaks of HCoV-229E, -NL63, and -HKU1 were detectable although they were obviously lower than the melting peak of HCoV-OC43 (Figure S2A). In particular, the peaks of HCoV-229E, -NL63, and -HKU1 increased along with the increase of primer concentration. When the concentrations of primers for HCoV-229E, -NL63, and -HKU1 were increased to 0.8  $\mu$ M, the melting peaks for four HCoVs appeared to have similar peak heights (Figure S2B). Therefore, 0.8, 0.8, 0.8 and 0.2  $\mu$ M of HCoV-229E, -NL63, -HKU1, and -OC43 specific primers were selected for the assay, respectively.



**Figure S1.** Comparison of the amplification effects between the original OC43 primers and the GC-modified OC43 primers. PC: positive control (OC43 standard strain VR-1558); NC: negative control; M: Marker.



**Figure S2.** Optimization of the primer concentrations for the multiple RT-qPCR assay. (A) Optimization of the multiple RT-qPCR assay using various concentrations of GC-modified HCoV-OC43 and other three HCoVs' primers; (B) Melting peaks of four HCoVs using the multiple RT-qPCR assay under optimized primer concentrations. PDs: primer dimers.





Figure S3. Standard curves of four HCoVs using the multiple RT-qPCR assay.





**Figure S5.** Melting curves of 9 HCoV-OC43, 2 HCoV-229E, and 9 HCoV-NL63 positive samples subjected to sequencing.