

Supplementary Table S1. Panels of SSR loci used.

Panel name	Locus	Primer name and final concentration used (μM)	Labelling
Multiplex1	MA4	MA4F (3)	PET
		MA4R (3)	
	MA7	MA7F (5)	FAM
		MA7R (5)	
	MA2-1	MA21F (5)	NED
		MA21R (5)	
Multiplex2	MA5-1	MA51F (5)	FAM
		MA51R (5)	
	MA14	MA14F (5)	FAM
		MA14R (5)	
	MA12	MA12F (5)	VIC
		MA12R (5)	
	MA19	MA19F (5)	FAM
		MA19R (5)	
	MA13	MA13F (5)	VIC
		MA13R (5)	