

Table S1. Oligonucleotide pairs used for the rhabdovirus detection.

Virus	Name	Sequence	Temperature
Chimay rhabdovirus	Chimay_L_F	5'-TCTTGAGCTTTGGTGGTCCG-3'	55°C
	Chimay_L_R	5'-CATCTTGTACACCCGCCCTT-3'	
Norway mononegavirus 1	Norway_L_F2	5'-ACACACTCCATCCTTTTGT-3'	50°C
	Norway_L_R2	5'-ACTACTCAAGCATCGGAATC-3'	
IRE/CTVM19-associated rhabdovirus	Rhabdo_Miass_L_1F	5'-GGGTTTGTGGTTAATTTGTC-3'	50°C
	Rhabdo_Miass_L_1R	5'-AGTGAGGACTGGATAAAAGA-3'	

Figure S1. Passage history of the IRE/CTVM19-associated rhabdovirus strains discussed in the study. Strain used in the current work marked in green.

