

Supplementary Materials

Pharmaceutics

Title:

The Effect of Cholesterol Content on the Adjuvant Activity of Nucleic-Acid-Free Lipid Nanoparticles

Authors:

Jessica Anindita^{1,2,†}, Hiroki Tanaka^{1,3,*,†}, Takuma Yamakawa², Yuka Sato², Chika Matsumoto¹, Kota Ishizaki², Taiji Oyama⁴, Satoko Suzuki⁵, Keisuke Ueda⁶, Kenjiro Higashi⁶, Kunikazu Moribe⁶, Kasumi Sasaki⁷, Yumika Ogura⁷, Etsuo Yonemochi⁷, Yu Sakurai¹, Hiroto Hatakeyama², and Hidetaka Akita^{1,3,*}

Affiliations

Jessica Anindita 1,2,†, **Hiroki Tanaka** 1,3,*,†, **Takuma Yamakawa** 2, **Yuka Sato** 2, **Chika Matsumoto** 1, **Kota Ishizaki** 2, **Taiji Oyama** 4, **Satoko Suzuki** 5, **Keisuke Ueda** 6, **Kenjiro Higashi** 6, **Kunikazu Moribe** 6, **Kasumi Sasaki** 7, **Yumika Ogura** 7, **Etsuo Yonemochi** 7, **Yu Sakurai** 1, **Hiroto Hatakeyama** 2 and **Hidetaka Akita** 1,3,*

¹ Laboratory of DDS Design and Drug Disposition, Graduate School of Pharmaceutical Sciences, Tohoku University, 6-3 Aoba Aramaki, Aoba-ku, Sendai City 980-8578, Miyagi, Japan

² Laboratory of DDS Design and Drug Disposition, Graduate School of Pharmaceutical Sciences, Chiba University, 1-8-1 Inohana, Chuo-ku, Chiba City 260-0856, Chiba, Japan

³ Center for Advanced Modalities and DDS, Osaka University, Suita 565-0871, Osaka, Japan

⁴ Sales Division, JASCO Corporation, 2967-5 Ishikawa, Hachioji City 192-8537, Tokyo, Japan; taiji.oyama@jasco.co.jp

⁵ Applicative Solution Lab Division, JASCO Corporation, 2967-5 Ishikawa, Hachioji City 192-8537, Tokyo, Japan

⁶ Laboratory of Pharmaceutical Technology, Graduate School of Pharmaceutical Sciences, Chiba University, 1-8-1 Inohana, Chuo-ku, Chiba City 260-0856, Chiba, Japan; keisuke@chiba-u.jp (K.U.)

⁷ Department of Physical Chemistry, School of Pharmacy and Pharmaceutical Sciences, Hoshi University, 2-4-41 Ebara, Shinagawa-ku, Shinagawa City 142-8501, Tokyo, Japan

* Correspondence: hiroki.tanaka.e1@tohoku.ac.jp (H.T.); hidetaka.akita.a4@tohoku.ac.jp (H.A.);

Tel.: +81-22-795-6831 (H.T. & H.A.)

† These authors contributed equally to this work.

Contents

S1. Detailed list of materials used in the study.

S2. Methods for supplementary figures

S2-1. Circular dichroism measurement

S2-2. Chemokine/cytokine measurement

S3. Supplementary figures

Figure S1. Correlation of lipid content and antibody production

Figure S2. Antibody production after separate administration of LNPs and antigens

Figure S3. Additional images of particle morphology observed using Cryo-TEM

Figure S4. The SAXS profile of Liposomes_A-D

Figure S5. Circular dichroism (CD) spectra

Figure S6. Dose-dependent chemokine/cytokine production

Figure S7. Chemokine/cytokine production of different nanoparticles

Supplementary Table S1. Detailed list of materials used in the study

Reagent	Size	Manufacturer	Product number
SM-102	100 mg	Cayman Chemical	33474
COATSOME MC-8080 (DSPC; 1,2-Distearoyl-sn-glycero-3-phosphocholine)	1 g	NOF CORPORATION	MC-8080
Cholesterol	5 g	Sigma Aldrich	C8667-5G
DMG-PEG (1,2-dimyristoyl-rac-glycero-3-methoxypolyethylene glycol) 2000	1 g	NOF CORPORATION	GM-020
Sodium Acetate (CH ₃ COONa)	500 g	nacalai tesque	C31118-75
MES (2-(<i>N</i> -morpholino)ethanesulfonic acid)	100 g	nacalai tesque	02442-44
HEPES (4-(2-hydroxyethyl)-1-piperazineethanesulfonic acid)	500 g	Dojindo	342-01375
D-PBS (-)	500 mL	nacalai tesque	14249-24
99.5% Ethanol	500 mL	nacalai tesque	14713-95
UltraPure™ DNase/RNase-Free Distilled Water (DDW)	500 mL	Invitrogen	10977-023
Albumin from chicken white egg	250 mg	Sigma Aldrich	A2512-250MG
Polyoxyethylene Sorbitan Monolaurate (Tween 20)	500 g	nacalai Tesque	35624-15
Fetal Bovine Serum (FBS)	500 mL	Gibco	10270
Goat anti-Mouse IgG-Fc Fragment Antibody HRP Conjugated	1 mL	Bethyl Laboratories	A90-131P
TMB Solution	100 mL	Merck Millipore	CL07-100ML
1mol/l Sulfuric Acid (H ₂ SO ₄)	500 mL	nacalai Tesque	95626-06
Sodium Hydrogen Carbonate (NaHCO ₃)	500 g	nacalai Tesque	31212-25
Sodium Carbonate (Na ₂ CO ₃)	500 g	nacalai Tesque	31310-35
LEGENDplex™ Mouse Anti-Virus Response Panel (13-plex) with V-bottom plate	100 tests	BioLegend	740622
Isoflurane Inhalation Anesthetic Solution	250 mL	Viatis	1119701G1092
AmiconUltra-4-100K Centrifugal Units	4 mL volume	Merck	UFC810096
AmiconUltra-4-100K Centrifugal Units	15 mL volume	Merck	UFC910096

S2. Methods for supplementary figures

S2-1. Circular dichroism (CD) measurement

LNP_A-D were prepared as described in the main text. The final concentration of samples was 6.6 mM of total lipid. Circular dichroism of these samples was measured using a J-1500 Circular Dichroism Spectrophotometer (JASCO Corporation, Japan). Parameters were as follows: measured range = 192-280 nm, scanning ratio = 50 nm/min, response = 4 sec, band width = 1 nm, accumulations = 4 times, and repeat = 3 times. The obtained CD spectrum was normalized by absorbance at 204 nm.

S2-2. Chemokine/cytokine measurement

Dispersed cholesterol particles, Liposome_C and LNP_A-D, were administered to BALB/c mice via s.c. injection (back of neck). Six hours after the injection, the blood was collected from the vena cava inferior using a 1 mL syringe and a 26G needle, and it then was moved to a 1.5 mL Eppendorf tube. The blood serum was obtained by incubating the blood at room temperature (RT) for 2 hours, followed by centrifugation (4 °C, 2,000G, 10 min). The chemokine/cytokine quantification was done using a LEGENDplex™ Mouse Anti-Virus Response Panel (13-plex) with a V-bottom Plate kit (#740622, BioLegend, California, US), following the manufacturer's instructions. These data were collected by NovoCyte Flow Cytometer Systems and analyzed using the Data Analysis Software Suite for LEGENDplex™. The chemokine/cytokine levels were determined from the mean fluorescence intensity (MFI) values.

Supplementary Figure S1




Figure S1. Correlation of lipid content and antibody production: Antibody production 14 days after boost vaccination correlated with **(a)** ionizable lipid SM-102, **(b)** DSPC, and **(c)** cholesterol content within the LNP. No significant correlation was found between antibody production and lipid content.

Supplementary Figure S2






Figure S2. Antibody production after separate administration of LNPs and antigens: The antibody production of LNP_A-D administered separately with recombinant OVA proteins is shown. Doses of 200 nmol of empty-LNPs were administered subcutaneously (s.c.) into BALB/c mice, followed by the administration (s.c.) of 10 μ g of rOVA 6 hours later. The treatment was performed twice with a 14-day interval. Later, the blood of each mouse was collected, and the serum was used for anti-OVA total IgG quantification by enzyme-linked immunosorbent assay (ELISA). The antibody production of LNPs did not depend on their ionizable lipid content. Mean \pm SD ($n = 1-5$).

Supplementary Figure S3


Dispersed cholesterol (98.5% chol)




Liposome_C (38.5% chol)




LNP_A (0% chol)



LNP_B (18.5% chol)



LNP_C (38.5% chol)



LNP_D (58.5% chol)




Figure S3. Additional images of particle morphology observed using Cryo-TEM: Additional Cryo-TEM images of dispersed cholesterol, Liposome_C, LNP_A, LNP_B, LNP_C, and LNP_D are shown. Each bar indicates 100 nm.

Supplementary Figure S4




Figure S4. The SAXS profile of Liposomes_A-D: The SAXS profiles of Liposome_A-D are shown. The SAXS profile of dispersed cholesterol particle is included for comparison.

Supplementary Figure S5




Figure S5. Circular dichroism (CD) spectra: (a) The CD spectra of LNP_A-D normalized at an absorbance wavelength of 204 nm; (b) The CD of LNP_A-C show a linear relationship with cholesterol content, while that of LNP_D shows a large absolute value compared with that of the other LNPs.

Supplementary Figure S6




Figure S6. Dose dependent chemokine/cytokine production: The chemokines/cytokines induced by LNP_C are shown with doses of 33.3 nmol, 66.6 nmol, and 100 nmol of total lipid concentration. The respective doses (33.3, 66.6, 100 nmol) of empty-LNPs were administered subcutaneously to the BALB/c mice. After 6 hours, the blood serum was collected and used to quantify the chemokine/cytokine levels. Mean with SD ($n = 3$), * $p < 0.05$, ** $p < 0.01$ (one-way ANOVA followed by Bonferroni against PBS).

Supplementary Figure S7




Figure S7. Chemokines/cytokines production of different nanoparticles: The chemokines/cytokines induced by dispersed cholesterol, Liposome_C, and LNP_A-D are shown. Doses of 100 nmol of empty LNPs were administered subcutaneously to the BALB/c mice; after 6 hours, the blood serum was collected and used to quantify the chemokine/cytokine levels. Mean with SD ($n = 3$), * $p < 0.05$, ** $p < 0.01$ (one-way ANOVA followed by Bonferroni against PBS).