

Table 1. List of antibodies used for immunocytochemistry and western immunoblotting,.

Protein	Catalog Numer, Company	Analysis, Dilution
E-cadherin	610181, BD Biosciences (Bath, UK)	WB, 1:1000, confocal microscopy 1:100,
β -actin	ab8227, Abcam (Cambridge, UK)	WB,1:2000
Integrin β 1	Ab1952P, Merck Millipore (Molsheim, France)	WB, 1:1000
β -catenin	610153, BD Biosciences (Bath, UK)	WB, 1:1000
Snail	3895, Cell Signaling Technology (Danvers, USA)	WB,1:1000
Alix	2171, Cell Signaling Technology (Danvers, USA)	WB,1:1000
Annexin V	8555, Cell Signaling Technology (Danvers, USA)	WB,1:1000
HSP70	4876, Cell Signaling Technology (Danvers, USA)	WB,1:1000
Cytochrome c	ab13575, Abcam (Cambridge, UK)	WB,1:1000
GAPDH	ab181602, Abcam (Cambridge, UK)	WB,1:20000
Glypican 1	16700-1-AP, proteintech (Manchester, UK)	WB, 1:600