

Table S1: Primers used in this study.

Target	Forward	Reverse
SARS-CoV-2 lab	TAGACATCATGCTAATGAGTACAGAT	TGAAGTCTTGTAAGTGTTCAG
	Probe ROX-GCTTATAACATGATGATCTCAGCTGGC-BHQ1	
HAdV-5 hexon	CCCITCGATGMTGCCCC	CIACGGGIACIAARCGCA
	Probe ROX-CCTGGCCCGAGATGTGCAT-BHQ1	
IAV M	The IAV primers and probe sequences are published as part of the AmpliSens® Influenza virus A/B-FL detection kit (Moscow, Russia, Reg.# RU #FSR 2009/05010).	