

Table S1: List of primers used in the Loop-mediated isothermal amplification to detect *invA* gene of *Salmonella*

PRIMER	Sequence (5'-3')
FIP	GACGACTGGTACTGATCGATAGTTTTTCAACGTTTCCTGCGG
BIP	CCGGTGAAATTATCGCCACACAAAACCCACCGCCAGG
F3	GGCGATATTGGTGTTCATGGGG
B3	AACGATAAACTGGACCACGG
LOOP F	GACGAAAGAGCGTGGTAATTAAC
LOOP B	GGGCAATTCGTTATTGGCGATAG

Primers developed by Hara-kudo (2005).

Hara-kudo, Y. Loop-Mediated Isothermal Amplification for the Rapid Detection of Salmonella. *FEMS Microbiol. Lett.*, **2005**, 253, 155–161. <https://doi.org/10.1016/j.femsle.2005.09.032>.