

Supplementary Information

Development of Robust and Standardized Cantilever Sensors Based on Biotin/Neutravidin Coupling for Antibody Detection. *Sensors* 2013, 13, 5273-5285

Jiayun Zhang ^{1,*}, Hans Peter Lang ¹, Felice Battiston ², Natalija Backmann ¹, Francois Huber ¹ and Christoph Gerber ¹

¹ Swiss Nano Institute, University of Basel, Klingelbergstrasse 82, 4056 Basel, Switzerland; E-Mails: hans-peter.lang@unibas.ch (H.P.L.); natalija.backmann@unibas.ch (N.B.); francois.huber@unibas.ch (F.H.); christoph.gerber@unibas.ch (C.G.)

² Concentris GmbH, Basel, Davidsbodenstrasse 63, 4056 Basel, Switzerland; E-Mail: battiston@concentris.ch

* Author to whom correspondence should be addressed; E-Mail: jiayun.zhang@unibas.ch; Tel.: +41-61-267-3769; Fax: +41-61-267-3784.

Figure A. Comparison of responses of IgG with multilayer which was formed using different functionalization methods. Orange bars: subsequent functionalization in the cantilever device by injecting 1.5 μ M NeutrAvidin solution into the measurement cell; violet bars: subsequent functionalization in 100 μ L of container by immersing the whole array in the same concentration of NeutrAvidin solution.

