Biosensors for the detection of interaction between Legionella pneumophila collagen-like protein and glycosaminoglycans

Han Su a, Shaopei Li a, Mauricio Terebiznik b, Cyril Guyard c and Kagan Kerman a,*

- ^a Department of Physical and Environmental Sciences, University of Toronto Scarborough, 1265 Military Trail, Toronto, M1C 1A4, ON, Canada.
- Department of Biological Sciences and Cell and Systems Biology, University of Toronto Scarborough, 1265
 Military Trail, Toronto, M1C 1A4, ON, Canada. terebiznik@utsc.utoronto.ca
- ^c BIOASTER Microbiology Technology Institute, 40 Avenue Tony Garnier, 69007 Lyon, France. cyril.guyard@bioaster.org
- * Corresponding author: kagan.kerman@utoronto.ca; Tel.: +1 416 287 7249

Supplementary Materials:





Figure S1: SPR responses of fucoidan and chondroitin Sulphate A on non-immobilized, Ni²⁺ activated NTA sensorchip. This data demonstrate that the non-specific interaction of the analytes with the Ni²⁺ activated NTA immobilized surface is negligible in the absence of His-rLcl.



 $Figure~S2: SPR~responses~of~Ni^{2+}~activation, followed~by~immobilization~of~His-rLcl.~HEPES~buffer~solution~was~injected~after~each~injection~to~remove~non-specifically~adsorbed~biomolecules.\\$